


WBC Industries Materials & Footcare  
Supplemental Catalog



LEATHER	3	DIGITAL CAPS™	17
NEOLON	4	BUNIONCARE GEL SLEEVE™	18
MICROCEL PUFF®	4	GEL BUNION™	18
WBC EVA FOAM SHEETS	5	GEL FOOT COVER™	18
QUICK-FORMABLES®	6	UNIVERSAL GEL STRAP™	19
PLASTAZOTE®	7	BALL OF FOOT GEL CUSHION™	19
SOFT SHEAR GEL SHEETING™	8	GEL METATARSAL PAD™	19
PORON® MEDICAL URETHANE SHEETS	8	ACHILLES HEEL PAD™	20
DUAL-LAMINATES	9	MALLEOLAR SLEEVE™	20
TRIPLE LAMINATES	9	SOFT SKIN HEEL SLEEVE™	20
ORTHO™CORK	10	BOOT BUMPER	21
CUSHION CORK	10	HEEL/ELBOW SLIP OVER™	21
THERMOCORK®	10	ARTHRITIC/DIABETIC GEL SOCK™	22
THERMOCORK® LITE	11	SOFTSOCK™	22
BIRKO-CORK	11	WONDERSPUR™	23
WBC BIO-KORK	11	WONDERCUP™	23
CARBOPLAST®	12	WONDERSPORT™	23
THERAPEUTIC FOOTCARE PRODUCTS	14	WONDERSOLE™	24
DIGITAL PADS ON A STRIP™	15	HEEL SPUR WITH REMOVABLE PLUG	24
DIGITAL PAD™	15	LOW PROFILE HEEL CUSHION	24
DIGISTRIP™	15	¾ LENGTH INSOLE	25
GEL TUBING™	16	FULL LENGTH INSOLE	25
GEL TOE CREST	16	FULL LENGTH INSOLE WITH TOE CREST	25
ALL GEL DIGITAL CAP	16	LOW PROFILE FABRIC-COVERED INSOLE	26
TOE SEPARATORS™	17	GEL METATARSAL PAD WITH TOE	26
TOE SPREADERS™	17	SPREADER™	

From early beginnings in the 1960s as the Westfield Brace Company, WBC Industries has evolved into one of the industry leaders in providing fabrication materials to the Orthotic and Prosthetic market.

Over the past few years, WBC has continued to expand its focus. In addition to supplying the O & P and Pedorthic markets with our own straps, buckles, hardware, and hook & loop, we've added an exciting new platform to our business: distribution! We now carry lines such as Silipos®, Rx Textiles, SPT Technology, & Waterhouse Leather, and we've expanded our selection of raw materials. We are now able to provide you with Plastazote, EVA, Poron, PPT, Corks, and Dual & Tri-Laminates from Acor and Apex.



Throughout this catalog you'll see many products for your everyday practice. With these new product offerings, we feel we are closer to meeting our goal: better serving you by satisfying your material and foot care product needs. We've never lost sight of what has always set WBC apart from other materials providers—our dedication to personal attention and focus on customer satisfaction.

CALL TODAY FOR A FREE LEATHER SWATCH SET



## LEATHER

- Sold by the hide.
- Average hide sizes:  
 Cow: 18-22 square feet  
 English Kips: 6-8 square feet  
 Molding & Russet: 8-12 square feet

Orthotic Materials

PRODUCT #	DESCRIPTION	PRODUCT #	DESCRIPTION
801	EXTRA LIGHT CREAM COW	148	BLACK SABLE
802	LIGHT CREAM COW	172	PEARL ELK (7/8 OZ.)
803	MEDIUM CREAM COW	181	LIGHT SMOKED ELK
901	EXTRA LIGHT PEARL COW	182	MEDIUM SMOKED ELK
902	LIGHT PEAL COW	183	HEAVY SMOKED ELK
903	MEDIUM PEARL COW	201	LIGHT BROWN ELK
297	BLACK COW	202	MEDIUM BROWN ELK
198	MEDIUM BLACK COW	203	HEAVY BROWN ELK
197	BROWN COW	1822	MEDIUM GREEN ELK
130	TAN LINING	1833	HEAVY GREEN ELK
195	WHITE COW	103	HEAVY BLACK ELK
300	BLUE COW	102	MEDIUM BLACK ELK
299	MAROON COW	101	LIGHT BLACK ELK
306	BEIGE CABRETTA	193	HEAVY WHITE ELK
307	CREAM CABRETTA	192	MEDIUM WHITE ELK
305	DARK BROWN CABRETTA	191	LIGHT WHITE ELK
1021	BLACK CALF	12	ENGLISH KIPS
1022	TAN CALF	154	MOLDING 5/6 OZ.
1023	BLACK HORSE	157	MOLDING 8/9 OZ.
1015	TOBACCO	704	RUSSET 5/6 OZ.
6241	TAN SUEDE	707	RUSSET 8/9 OZ.
161	BLACK LUXOR	804	HEAVY CREAM COW
163	NAVY LUXOR	1013	PEARL APRON SPLITS
164	BONE LUXOR	128	BLACK SADDLE SIDE
162	BROWN LUXOR	127	BROWN SADDLE SIDE
191A	WHITE ATHLETIC	129	CEDAR SADDLE SIDE
147	IVORY SABLE	146	HONEY SABLE

## NEOLON

New & improved firmer neoprene – for more active patients. Colored nylon top cover with abraded bottom side.

PROD. #	THICKNESS	COVER COLOR	SIZE
416	1/16"	ROYAL BLUE	41" X 51"
426	1/16"	ROYAL BLUE	82" X 51"
418	1/8"	ROYAL BLUE	41" X 51"
428	1/8"	ROYAL BLUE	82" X 51"
516	1/16"	BLACK	41" X 51"
526	1/16"	BLACK	82" X 51"
518	1/8"	BLACK	41" X 51"
528	1/8"	BLACK	82" X 51"
616	1/16"	WINE	41" X 51"
626	1/16"	WINE	82" X 51"
618	1/8"	WINE	41" X 51"
628	1/8"	WINE	82" X 51"
718	1/8"	GREEN	41" X 51"
728	1/8"	GREEN	82" X 51"




## MICROCEL PUFF®

- 35 Durometer EVA foam
- Excellent shock absorber – resists bottoming out
- Heat moldable with excellent shape retention
- Dense, durable and grinds easily with a clean finish
- Excellent for firm support in arch supports
- Molding Temperature: 275-300° F - Heat Time: 2-3 Minutes
- Sheet size: 37 x 60"
- Thicknesses: 1/16" 1/8" 3/16" 1/4" 3/8" 1/2" 5/8" 3/4" 1"

STYLES	COLORS
SOLIDS:	BLACK, BLUE, FOREST GREEN, GREEN, PINK, PURPLE, WHITE, YELLOW
CONFETTI:	BLUE / BLACK / AQUA YELLOW / RED / BLUE / GREEN
SWIRL:	GREEN / PURPLE / WHITE PINK / PURPLE / WHITE

*25, 55, & 65 Durometer Microcel Puff are also available.*





Purple

Red

Watercolor

Blue

Flesh

## WBC EVA FOAM SHEETS

Used for foot orthoses and prosthetic devices, orthotic posting and heel lifts. Can be used as a base, middle layer, or top cover for many foot orthoses. Highly resilient and washable.

Orthotic Materials

PRODUCT #	COLOR	THICKNESS	DUROMETER	SHEET SIZE
169-10-02	PURPLE SWIRL	2.0 MM	35	25" X 43"
169-10-03	PURPLE SWIRL	3.0 MM	35	25" X 43"
169-10-04	PURPLE SWIRL	4.0 MM	35	25" X 43"
169-20-02	RED	2.0 MM	35	25" X 43"
169-20-03	RED	3.0 MM	35	25" X 43"
169-20-04	RED	4.0 MM	35	25" X 43"
169-30-02	WATERCOLOR	2.0 MM	35	25" X 43"
169-30-03	WATERCOLOR	3.0 MM	35	25" X 43"
169-30-04	WATERCOLOR	4.0 MM	35	25" X 43"
169-40-01	WHITE	1.5 MM	30	39" X 39"
169-40-03	WHITE	3.0 MM	30	39" X 39"
169-40-04	WHITE	4.0 MM	30	39" X 39"
169-50-30	BLUE	3.0 MM	30	42" X 50"
169-50-04	BLUE	4.0 MM	30	39" X 39"
169-50-40	BLUE	4.0 MM	30	42" X 50"
169-60-30	FLESH	3.0 MM	30	42" X 50"
169-60-04	FLESH	4.0 MM	30	39" X 39"
169-70-01	BLUE	1.5 MM	40	39" X 39"
169-70-30	BLUE	3.0 MM	40	38" X 46"
169-70-40	BLUE	4.0 MM	40	40" X 46"

*Also available in 5mm, 6mm, 10mm thicknesses.*

## QUIKFORMABLES®

Designed to cut fabrication time in half, QUIKformables® help you produce truly custom-made orthoses that are indistinguishable from orthoses made by traditional methods. QUIKformables® aid in the cast modification process, making this chore simple, neat and fast.

- Custom orthoses in a fraction of the time
- Available in any material combination
- Just heat, vacuum form, and grind

### Footbed Material (choose one):

PRODUCT #	DESCRIPTION
QFAB-01	1/4" PUFF
QFAB-02	1/4" MULTICORK
QFAB-03	PINK PLASTAZOTE + BLUE PORON + PUFF
QFAB-04	PUFF + BLUE PORON + PUFF
QFAB-05	SLOW RECOVERY PORON + PUFF
QFAB-15	P-CELL + BLUE PORON + PUFF



### Posting Material:

ITEM #	DESCRIPTION
01	MULTICORK
02	PUFF
03	EFM

### Sizes Available:

<b>WOMEN</b>	4-5	6-7	8-9	10-11	12-13	14-15	16-17
<b>MEN</b>	2-3	4-5	6-7	8-9	10-11	12-13	14-15

*Also Available:*

### BLOCKERS from Acor and Apex

- Reduce your materials costs and valuable fabrication time with these pre-laminated, pre-skived blockers.
- Blockers are tapered, thinner at the toe and thicker at the heel, allowing for post-heating modifications such as grinding or posting.
- Blockers combine state-of-the-art orthotic materials that, when heated, allow you to custom fabricate custom foot orthoses quickly and easily.

## PLASTAZOTE®

- Closed cell, cross-linked polyethylene foam
- Heat formable at 275° F
- Comfortable, lightweight, conformable
- Great for diabetic & arthritic usage



### #1 MEDIUM: 15 SHORE A DUROMETER

COLOR:	THICKNESSES:	SHEET SIZE:
BLACK	1/16, 1/8, 3/16, 1/4	30 X 40"
WHITE	1/16, 1/8, 3/16, 1/4, 3/8, 1/2, 3/4, 1	30 X 40"
PINK	1/16, 1/8, 3/16, 1/4, 3/8, 1/2, 3/4, 1	30 X 40"
RED	1/8, 1/4, 1/2	30 X 40"
YELLOW	1/8, 1/4, 1/2	30 X 40"
BLUE	1/8, 1/4, 1/2	30 X 40"
GREEN	1/8, 1/4, 1/2	30 X 40"

### #2 FIRM: 30 SHORE A DUROMETER

COLOR:	THICKNESSES:	SHEET SIZE:
BLACK	1/16, 1/8, 3/16, 1/4, 3/8, 1/2, 3/4, 1	34 X 34"
WHITE	1/16, 1/8, 3/16, 1/4, 3/8, 1/2, 3/4, 1	34 X 34"

### #3 RIGID: 60 SHORE A DUROMETER

COLOR:	THICKNESSES:	SHEET SIZE:
BLACK	1/16, 1/8, 3/16, 1/4, 3/8, 1/2	25 X 35"
WHITE	1/16, 1/8, 3/16, 1/4, 3/8, 1/2	25 X 35"


### PLASTAZOTE® CASE QUANTITIES

1" = 3 SHEETS	3/4" = 4 SHEETS	1/2" = 6 SHEETS
3/8" = 8 SHEETS	1/4" = 12 SHEETS	3/16" = 16 SHEETS
1/8" = 24 SHEETS	1/16" = 24 SHEETS	

**SOFT SHEAR GEL SHEETING™**

- 3mm thick sheet of Silipos® gel, available either lined on both sides, or just one side
- Fabric cover chosen for stretch qualities, excellent pressure and shear relieving properties
- Fabric increases durability of gel & allows sheet to be glued
- One sheet per package

PRODUCT #	SIZE	DESCRIPTION
15605	8 X 8" (20 CM X 20 CM)	FABRIC 1 SIDE
15615	8 X 8" (20 CM X 20 CM)	FABRIC 2 SIDES
15705	12 X 16" (30.5 CM X 40.5 CM)	FABRIC 1 SIDE
15715	12 X 16" (30.5 CM X 40.5 CM)	FABRIC 2 SIDES
15805	16 X 20" (40.5 CM X 51 CM)	FABRIC 1 SIDE
15815	16 X 20" (40.5 CM X 51 CM)	FABRIC 2 SIDES

**PORON® MEDICAL URETHANE SHEETS**

Open-cell urethane is breathable, odorless, washable and non-sensitizing. Will not bottom out over time, and easily bonds to other materials.

STYLES:	SHEET SIZES:	COLORS:	THICKNESSES:
SMOOTH ON BOTH SIDES	12 X 42" 18 X 42" 36 X 42"	BLUE TAN	1/16, 1/8, 1/4, 3/8, 1/2 1/8, 3/8
ABRADED ON ONE SIDE	12 X 42" 18 X 42" 36 X 42"	BLUE TAN	1/16, 1/8, 1/4, 3/8 1/16, 1/8, 1/4, 3/8
ABRADED ON BOTH SIDES	12 X 42" 18 X 42" 36 X 42"	BLUE TAN	1/16, 1/8, 1/4, 3/8 1/16, 1/8, 1/4,
ABRADED ON ONE SIDE, NAVY CLOTH ON THE OTHER	12 X 42" 18 X 42" 36 X 42"	BLUE	1/16, 1/8
SCRIM/ABRADED	12 X 42" 18 X 42" 36 X 42"	BLUE	1/16, 1/8

## DUAL-LAMINATES:

### PLASTAZOTE® + PORON®

PRODUCT #	DESCRIPTION	THICK	SHEET SIZE
PZPN-00014	1/8" #1 PINK PLASTAZOTE + 1/16" BLUE PORON	3/16"	37" X 14"
PZPN-00006	1/8" #1 PINK PLASTAZOTE + 1/16" BLUE PORON	3/16"	37" X 28"
PZPN-00015	1/8" #1 PINK PLASTAZOTE + 1/8" BLUE PORON	1/4"	37" X 14"
PZPN-00007	1/8" #1 PINK PLASTAZOTE + 1/8" BLUE PORON	1/4"	37" X 28"
PZPN-00013	1/4" #1 PINK PLASTAZOTE + 1/8" BLUE PORON	3/8"	37" X 14"
PZPN-00005	1/4" #1 PINK PLASTAZOTE + 1/8" BLUE PORON	3/8"	37" X 28"

Also Available:

P-CELL™ + PORON®

P-CELL™ + SLOW RECOVERY PORON® (SRP)

PLASTAZOTE® + SLOW RECOVERY PORON® (SRP)

ORTHACOR™

## TRIPLE LAMINATES:

### P-CELL™ + PORON® + MICROCEL PUFF®

PRODUCT #	DESCRIPTION	THICK	SHEET SIZE
MDUR-00083	1/8" P-CELL + 1/16" BLUE PORON + 1/8" BLUE PUFF	5/16"	37" X 30"
MDUR-00084	1/8" P-CELL + 1/16" BLUE PORON + 1/4" BLUE PUFF	7/16"	37" X 30"
MDUR-00085	1/4" P-CELL + 1/16" BLUE PORON + 1/8" BLUE PUFF	7/16"	37" X 30"
MDUR-00086	1/8" P-CELL + 1/8" BLUE PORON + 1/8" BLUE PUFF	3/8"	37" X 30"

### PLASTAZOTE® + PORON® + MICROCEL PUFF®

PRODUCT #	DESCRIPTION	THICK	SHEET SIZE
MDUR-00001	1/8" #1 PINK PLASTAZOTE + 1/16" BLUE PORON + 1/8" BLUE PUFF	5/16"	37" X 30"
MDUR-00002	1/8" #1 PINK PLASTAZOTE + 1/16" BLUE PORON + 1/4" BLUE PUFF	7/16"	37" X 30"
MDUR-00003	1/4" #1 PINK PLASTAZOTE + 1/16" BLUE PORON + 1/8" BLUE PUFF	7/16"	37" X 30"
MDUR-00058	1/8" #1 PINK PLASTAZOTE + 1/8" BLUE PORON + 1/8" BLUE PUFF	3/8"	37" X 30"

## TRI-LAMINATES

PRODUCT #	DESCRIPTION	THICK	SHEET SIZE
TRLM-00001	1/16" PINK PUFF + 1/16" BLUE PORON + 1/8" PINK PUFF	1/4"	37" X 30"
TRLM-00007	1/16" BBA PUFF + 1/16" BLUE PORON + 1/8" PINK PUFF	1/4"	37" X 30"
TRLM-00003	1/16" BLUE PUFF + 1/16" BLUE PORON + 1/8" BLUE PUFF	1/4"	37" X 30"
TRLM-00005	1/16" WHITE PUFF + 1/16" BLUE PORON + 1/8" WHITE PUFF	1/4"	37" X 30"

Also Available:

P-CELL™ + SLOW RECOVERY PORON® (SRP) + MICROCEL PUFF®

PLASTAZOTE® + SLOW RECOVERY PORON® (SRP) + MICROCEL PUFF®

TRILAM + SLOW RECOVERY PORON® (SRP)

MICROCEL PUFF® + EFM + MICROCEL PUFF®

**ORTHO™CORK**

A cork sheet designed to be light-weight and strong, which can be shaped, ground, and trimmed. Great for elevations, shoe inserts, and orthotic filler.

PRODUCT #	DESCRIPTION	THICKNESS	SHEET SIZE
8011-1	ORTHOCORK	1/16"	18" X 36"
8011-2	ORTHOCORK	1/8"	18" X 36"
8011-3	ORTHOCORK	3/16"	18" X 36"
8011-4	ORTHOCORK	1/4"	18" X 36"
8011-5	ORTHOCORK	3/8"	18" X 36"
8011-6	ORTHOCORK	1/2"	18" X 36"

**CUSHION CORK**

A cork-rubber mix that provides flexibility and durability. Easy to grind – use for a durable base or a reinforcement for foot orthotics.

PRODUCT #	DESCRIPTION	THICKNESS	SHEET SIZE
8012-1	CUSHION CORK	1/16"	18" X 36"
8012-2	CUSHION CORK	1/8"	18" X 36"
8012-3	CUSHION CORK	3/16"	18" X 36"
8012-4	CUSHION CORK	1/4"	18" X 36"
8012-5	CUSHION CORK	3/8"	18" X 36"
8012-6	CUSHION CORK	1/2"	18" X 36"

**THERMOCORK®**

- A thermoplastic orthotic cork that is heat-formable to 170° F
- Can be formed directly to the positive mold
- Strong & flexible, it will not pack or bottom out
- Cover with Plastazote, STS, or leather to customize your orthoses to the appropriate medical condition
- Easy to grind & trim.

PRODUCT #	DESCRIPTION	THICKNESS	SHEET SIZE
8013-2	THERMOKORK	1/8"	36" X 44"
8013-3	THERMOKORK	3/16"	36" X 44"
8013-4	THERMOKORK	1/4"	36" X 44"
8013-5	THERMOKORK	3/8"	36" X 44"

**THERMOCORK® LITE**

- A lightweight thermoplastic blend of cork and EVA
- Provides greater flexibility & cushioning
- Easy to mold & grind
- 55 Durometer

PRODUCT #	DESCRIPTION	THICKNESS	SIZE
8014A-1	THERMOCORK LITE	1/16"	18" X 36"
8014A-2	THERMOCORK LITE	1/8"	18" X 36"
8014A-3	THERMOCORK LITE	3/16"	18" X 36"
8014A-4	THERMOCORK LITE	1/4"	18" X 36"
8014A-5	THERMOCORK LITE	3/8"	18" X 36"

**BIRKO-CORK**

Made from natural cork granules mixed with nylon, Birko-Cork is porous, lightweight, flexible, and heat formable. Material becomes pliable but not sticky at 174°F. At 212°F it is slightly sticky and good for molding. At 250°F it is sticky, excellent for molding, and adheres well to other materials. Leftovers from each job may be kneaded together for future use.

PRODUCT #	DESCRIPTION	THICKNESS	SHEET SIZE	DENSITY
13105	BIRKO-CORK	3/16"(5.0 MM)	12" X 18"	SOFT
13010	BIRKO-CORK	3/8"(10.0 MM)	12" X 18"	SOFT
13305	BIRKO-CORK	3/16"(5.0 MM)	12" X 18"	HARD
13310	BIRKO-CORK	3/8"(10.0 MM)	12" X 18"	HARD

**WBC BIO-KORK**

- A traditional blend of EVA & cork materials
- 60-65 durometer, vacuum formable, and easy to grind
- Heat formable at 300°F

PRODUCT #	DESCRIPTION	THICKNESS	SHEET SIZE
33003	WBC BIO-KORK	1/8"(3.0 MM)	32" X 40"
33004	WBC BIO-KORK	5/32"(4.0 MM)	32" X 40"
33006	WBC BIO-KORK	1/4"(8.0 MM)	32" X 40"
33008	WBC BIO-KORK	5/16"(8.0 MM)	32" X 40"
33010	WBC BIO-KORK	3/8"(10.0 MM)	32" X 40"
33012	WBC BIO-KORK	1/2"(12.0 MM)	32" X 40"

## CARBOPLAST® SHEETS

Carboplast® is an easy to mold, light-weight, and durable material made of high-strength carbon and glass fibers in a tough, unique polypropylene matrix.

- Very light-weight, ultra thin
- Guaranteed not to crack for life
- Molds deep heel cups using conventional equipment
- 4 rigidities to meet all patients' demands
- Molding Temperature: 400° F

### STANDARD CARBOPLAST® MATERIAL

PROD. #	NAME	SIZE	THICKNESS
80271	CARBOPLAST - FLEXIBLE	A - F*	1.5MM
80273	CARBOPLAST - SEMI-RIGID	A - F*	2.0MM
80275	CARBOPLAST - RIGID	A - F*	2.5MM
80277	CARBOPLAST - X-RIGID	A - F*	2.9MM

### ICE CARBOPLAST® MATERIAL

PROD. #	NAME	SIZE	THICKNESS
80271I	CARBOPLAST - FLEXIBLE	A - F*	1.5MM
80273I	CARBOPLAST - SEMI-RIGID	A - F*	2.0MM
80275I	CARBOPLAST - RIGID	A - F*	2.5MM
80277I	CARBOPLAST - X-RIGID	A - F*	2.9MM

### FIRE CARBOPLAST® MATERIAL

PROD. #	NAME	SIZE	THICKNESS
80271F	CARBOPLAST - FLEXIBLE	A - F*	1.5MM
80273F	CARBOPLAST - SEMI-RIGID	A - F*	2.0MM
80275F	CARBOPLAST - RIGID	A - F*	2.5MM
80277F	CARBOPLAST - X-RIGID	A - F*	2.9MM

### EXTREME CARBOPLAST® MATERIAL

PROD. #	NAME	SIZE	THICKNESS
80271X	CARBOPLAST - FLEXIBLE	A - F*	1.5MM
80273X	CARBOPLAST - SEMI-RIGID	A - F*	2.0MM
80275X	CARBOPLAST - RIGID	A - F*	2.5MM
80277X	CARBOPLAST - X-RIGID	A - F*	2.9MM

\*Carboplast® Sheet Sizes = A - 3.5" x 7.5", B - 4" x 8", C - 5" x 9", D - 4" x 12", E - 18" x 36", F - 36" x 36"

Standard

Ice

Fire

Extreme



## CARBOPLAST® Thermoformable Plates

- Used to limit motion for conditions such as turf toe or hallux rigidus
- Easy to fit and finish – can be customized by grinding or heat molding
- Can be heated and reformed by using convection ovens at 400° F for 3 - 7 minutes
- Guaranteed not to crack or lose its shape or material properties

PRODUCT #	DESCRIPTION
4006 (1-5)	CARBOPLAST CONTOURED PLATES - SIZES: S,M,L,XL
4007 (1-5)	CARBOPLAST FLAT PLATES - SIZES: XS,S,M,L,XL
4008 (3-5)	CARBOPLAST CONTOURED MORTON'S TOE PLATES (L & R)
4009 (3-5)	CARBOPLAST FLAT MORTON'S TOE PLATES

**Thermoformable Plates:**  
Please specify **Flexible (1)**,  
**Semi-Rigid (3)**, or **Rigid (5)**



## CARBOPLAST® Orthotics

- Pre-fabricated foot orthotics that offer rigid support while remaining light-weight and extra thin
- Will not crack, even under severe pressure, and will fit comfortably in all footwear styles
- Available with and without metatarsal support
- Can be heated and remolded for a custom fit

PROD. #	DESCRIPTION
1445-0	CARBOPLAST SHAFFER W/MET – WOMEN 4 - 11, MEN 6 - 15
1445-1	CARBOPLAST SHAFFER NO MET - WOMEN 4 - 11, MEN 6 - 15



1445-0



1445-1

## CARBOPLAST® Sole Bars

- Light-weight, rigid and used when motion control is needed
- Will limit a shoe's flexion
- Can be heat adjusted

PRODUCT #	DESCRIPTION
40051	CARBOPLAST SOLE BAR - 2MM X 1" X 8", 9", 10", 11", 12"



## Therapeutic Footcare Products

Your patients' feet go through a lot. Why not give them the best in affordable foot care products? We're proud to offer the latest Silipos® Gel technology in our Therapeutic Foot Care line.

Silipos® Gel contains a medical-grade mineral oil which gradually diffuses into the skin to moisturize, soften and condition. The gel's fluid-like properties allow it to move with the skin and absorb shear forces, pressure, and friction. It's non-toxic, hypoallergenic, and antimicrobial - the perfect choice for your discerning patients.

### DIGITAL PADS ON A STRIP™

- Elasticized tubing is easy to cut using sharp scissors
- Once cut to size, the tube can be turned inside-out
- End result is virtually the same as Digital Pads™

PROD #	SIZE	USAGE	QUANTITY	SOLD
1105	MINI	5TH TOE	17 PADS/STRIP	2/PKG
1106	S/M	LESSER TOES OR FINGERS	16 PADS/STRIP	2/PKG
1107	L/XL	HALLUX OR THUMB	14 PADS/STRIP	2/PKG



### DIGITAL PAD™

- Protects toes from friction and footwear pressure
- Use where corns, calluses, or hammer toes cause discomfort

PRODUCT #	SIZE	USAGE	SOLD
10055	MINI	5TH TOE	12/PKG
10065	S/M	LESSER TOES OR FINGERS	12/PKG
10075	L/XL	HALLUX OR THUMB	12/PKG



### DIGISTRIP™

- Half-coated with Silipos® polymer gel
- Very cost effective
- Same applications as the digital pad but in a strip, allowing you to choose the length needed
- Approximately 24" (60cm) long

PRODUCT #	SIZE	SOLD
1100	NARROW	2 STRIPS/PKG
1101	MEDIUM	2 STRIPS/PKG
1102	WIDE	2 STRIPS/PKG



## GEL TUBING™

- Fully lined with Silipos® polymer gel
- All-around digital protection
- Uses less room than the conventional foam tubing
- Gel will not bottom out as foam materials do
- Can help reduce scar tissue
- Can be cut to size required for toes or fingers

PRODUCT #	STYLE	SIZE	SOLD
10865	RIBBED KNIT FABRIC	NARROW, 3/4" X 6"	2/PKG
10875	RIBBED KNIT FABRIC	WIDE, 1" X 6"	2/PKG
10855	MESH	NARROW, 3/4" X 6"	3/PKG
10845	MESH	WIDE, 3" X 10"	1/PKG



## GEL TOE CREST

- 100% medical grade silicone absorbs shock & vibration
- Provides even pressure distribution
- Helps reduce pressure and pain caused by hammer, claw, and mallet toes
- Will not bottom out over time

PRODUCT #	SIZE	STYLE	SOLD
10415	SMALL	LEFTS	3 IN A POLY BAG
10425	SMALL	RIGHTS	3 IN A POLY BAG
10435	LARGE	LEFTS	3 IN A POLY BAG
10445	LARGE	RIGHTS	3 IN A POLY BAG
10515	SMALL	LEFT	1 IN A BOX
10525	SMALL	RIGHT	1 IN A BOX
10535	LARGE	LEFT	1 IN A BOX
10545	LARGE	RIGHT	1 IN A BOX



## ALL GEL DIGITAL CAP

- Fully lined with Silipos® gel to moisturize, condition and soften the skin
- The all-gel digital cap is ideal for easing pressure and friction caused by corns, blisters, calluses, hammertoes, and in-grown nails
- Virtually transparent, they are ideal for open toe shoes/sandals
- Average cap length: 2.5" (6.5cm)
- One size fits all. Washable and reusable

PRODUCT #	SIZE	SOLD
10495	S/M	2/PKG
10505	M/L	2/PKG



## TOE SEPARATORS™

- Anatomically molded to fit comfortably between toes, reducing shearing and inter-digital pressure
- Can be trimmed if required
- Cost-effective alternative to silicone compounds

PRODUCT #	SIZE	SOLD
11605	SMALL	15/BAG
11615	LARGE	15/BAG



## TOE SPREADERS™

- Thicker than the Toe Separators™
- Due to the polymer's soft conforming properties, these spreaders are extremely comfortable
- Used as an inter-digital wedge for bunion deformity, overlapping toes, or toe drift
- Can be trimmed with scissors

PRODUCT #	SIZE	SOLD
11505	SMALL	4/BAG
11515	MEDIUM	4/BAG
11525	LARGE	4/BAG



## DIGITAL CAPS™

- Fully lined with Silipos® polymer gel, the digital cap is ideal for all-round digital protection
- Widely used in sports for pressure and friction protection
- Can be used to protect over/under-lapping toes
- Available in two materials—stretchable fabric offering a mild compression depending on sizes chosen, and mesh material that is more elasticized.
- Average cap length is 2.5" – fabric will not fray when cut

PRODUCT #	STYLE	SIZE	SOLD
10195	RIBBED KNIT FABRIC	MINI	6/PKG
10205	RIBBED KNIT FABRIC	S/M	6/PKG
10206	RIBBED KNIT FABRIC	L/XL	6/PKG
10235	MESH	ONE SIZE	8/PKG



### BUNIONCARE GEL SLEEVE™

- Comfortable sleeve cushions, pads, and protects
- Anatomically designed gel pad conforms to the shape of the bunion
- Alleviates friction and pressure
- Helps reduce scarring after surgery

PRODUCT #	SIZE	SOLD
10325	S/M	1 PAIR/PKG
10335	L/XL	1 PAIR/PKG



### GEL BUNION™

- Protects bunion joint using a concave gel pad
- Soft foam shell for additional pressure absorption
- Toe loop can be cut at the spoon end for a custom comfortable fit
- Gel pad comes separate from shell allowing exact positioning for each patient
- Can be worn with dress shoes

PRODUCT #	SIZE	SOLD
10085	S	5/PKG
10095	L	5/PKG



### GEL FOOT COVER™

- Made from a thin layer of gel, covered on both sides with fabric
- Offers complete toe protection
- Extends to under the metatarsal heads
- Can be worn comfortably in everyday footwear, even under hosiery
- Ideal if an activity places stress on toes

PROD. #	SIZE	SOLD
10275	SMALL - 3.75" X 4.9" (9.5 CM X 12.5 CM)	1 PAIR/PKG
10285	MEDIUM - 4.0" X 5.3" (10 CM X 13.5 CM)	1PAIR/PKG
10295	LARGE - 4.3" X 5.9" (11 CM X 15 CM)	1PAIR/PKG



### UNIVERSAL GEL STRAP™

- Surgical-grade elastic sleeve with a gel pad that covers all five metatarsal heads
- Helps reduce pain and discomfort
- Serves as a soft tissue supplement to fat pad atrophy
- Pad is available uncovered for maximum benefit of the mineral oil based gel
- Pad is also available covered to provide an extended life when callusing is not a problem
- By cutting the stitch, can also be used on the instep, wrist or hand



PRODUCT #	SIZE	SOLD
10155	S/M (GEL UNCOVERED)	1/PKG
10165	L/XL (GEL UNCOVERED)	1/PKG
10175	S/M (GEL COVERED)	1/PKG
10185	L/XL(GEL COVERED)	1/PKG

### BALL OF FOOT GEL CUSHION™

- Excellent for relieving the burning sensation under the metatarsal heads of the foot
- Protects against shear forces, abrasion, friction, and pressure
- May be used as a soft tissue supplement for Metatarsalgia
- Ideal for calluses and small neuromas
- Gel thickness is 2mm



PRODUCT #	SOLD
10025	2PCS./PKG

### GEL METATARSAL PAD

- 100% medical grade silicone absorbs shock, vibration, and provides even pressure distribution
- Anatomically designed for comfort and support
- Helps reduce the burning sensation under your metatarsal heads.
- Ideal for conditions such as Fat Pad Atrophy, Metatarsalgia, Morton's Neuroma, etc.

PRODUCT #	SIZE	SOLD
10465	ONE SIZE FITS ALL	1 PAIR/PKG



### ACHILLES HEEL PAD™

- A soft conforming gel pad protects the posterior of the heel from pressure and friction
- Provides gentle compression and dissipation of pressure around the inflamed area
- Effective way to protect the back of the heel from pressure and friction caused by sports boots: Ski, In-Line, Hockey & Ice Skating
- Reduces pain caused by footwear pressure, particularly on any irregular heel bumps

PRODUCT #	SIZE	SOLD
10385	S/M	1/PKG
10395	L/XL	1/PKG

GEL THICKNESS: 1/8" (3 MM)



### MALLEOLAR SLEEVE™

- Two large discs of polymer gel molded on both medial and lateral malleolar, reducing pressure and friction to the area
- Can be used pre or post injury
- Provides gentle compression to ankle bones
- Effective way to protect ankles from sports boots: Ski, In-Line, Hockey, & Ice Skating
- Reduces abrasion and shear forces

PRODUCT #	SIZE	SOLD
15005	ONE SIZE	1/PKG



### SOFT SKIN HEEL SLEEVE™

- Gel pad anatomically molded to the shape of the heel
- Will moisturize and soften dry, hard skin, helping to reduce heel fissures
- Occlusion has been shown to be beneficial for psoriasis and eczema
- Can be used in conjunction with skin moisturizing cream

PRODUCT #	SIZE	SOLD
15265	ONE SIZE	1/PKG



### BOOT BUMPER™

- Conforming polymer gel pads that protect the front and back of the ankle from pressure, pinching, and friction
- Helps protect the foot and ankle from abrasions caused by athletic footwear, without the added bulk
- Effective way to protect ankles from sports boots:  
Ski – In-Line – Hockey – Ice Skating

PRODUCT #	SIZE	SOLD
1770	S/M	1/PKG
1772	L/XL	1/PKG



### HEEL/ELBOW SLIP OVER™

- Helps the healing of dry, cracked heels or elbows
- Used to help prevent heel ulceration
- Gel pad can be placed on dorsal side of foot to absorb shoe lace bite
- Pad can be positioned anywhere around the foot/heel

PRODUCT #	SIZE	SOLD
15225	S/M	1/PKG
15235	L/XL	1/PKG



### ARTHRITIC/DIABETIC GEL SOCK™

- Single layer sock with 1/16" (2mm) of gel throughout the sole
- Helps prevent foot ulcerations
- Seamless toe
- Non-restrictive cuff
- Use with or without diabetic shoes
- Optimum protection for the neuropathic foot
- Helps reduce friction, abrasion, shear forces & callus formation
- Serves as a soft tissue supplement for fat pad atrophy
- One pair per pack



PRODUCT #		U.S.	OTHER		SHOE SIZE
WHITE	BLACK	SOCK SIZE	SOCK SIZE	CM	
1701	1706	7-9	37-40	25-27	FEMALE SIZE 6 OR UNDER
1702	1707	9-11	41-43	27-29	FEMALE OVER SIZE 6/ MALE UNDER SIZE 9
1703	1708	10-13	44-46	28-31	MALE SIZE 9 AND OVER

### SOFTSOCK™

- A double layer sock with a 1/4" (6mm) gel insole between the sole layers
- For any condition that requires a soft, friction-absorbing foot bed
- Indicated for neuropathic & rheumatoid arthritis sufferers
- Fluid-like nature of the gel massages the foot while walking
- Reduces skeletal shock on heel strike
- Serves as a soft tissue supplement for fat pad atrophy
- One pair per pack



PRODUCT #	COLOR	SIZE	SOCK SIZE	SOLD
10605	WHITE	SMALL	7-9	1 PAIR/PKG
10615	WHITE	MEDIUM	9-11	1 PAIR/PKG
10625	WHITE	LARGE	10-13	1 PAIR/PKG
10635	WHITE	X-LARGE	13-16	1 PAIR/PKG
10705	BLACK	SMALL	7-9	1 PAIR/PKG
10715	BLACK	MEDIUM	9-11	1 PAIR/PKG
10725	BLACK	LARGE	10-13	1 PAIR/PKG
10735	BLACK	X-LARGE	13-16	1 PAIR/PKG

## WONDERZORB™ Dual Density Orthotics

Maximum support and shock absorption using two densities of medical-grade silicone – clear silicone with a softer density colored silicone dots – provide additional relief to sensitive areas.

### WONDERSPUR™

- Heel Spurs
- Knee, back & ankle pain
- Plantar Fasciitis
- Fat Pad Atrophy

PRODUCT #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4500	SMALL	2-4	4-5
4501	MEDIUM	5-7	6-9
4502	LARGE	8-10	10-12
4503	X-LARGE	11 & UP	12 & UP



### WONDERCUP™

- Bursitis
- Achilles Tendonitis
- Pressure Ulcers
- Plantar Fasciitis

PRODUCT #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4504	SMALL	2-4	4-5
4505	MEDIUM	5-7	6-9
4506	LARGE	8-10	10-12
4507	X-LARGE	11 & UP	12 & UP



### WONDERSPORT™

- Metatarsalgia
- Fat Pad Atrophy
- Diabetic Feet
- Elongated Metatarsals
- Morton's Neuroma

PRODUCT #	MEN'S SIZES	PRODUCT #	WOMEN'S SIZES
4508	5-7	4513	4-5
4509	8-9	4514	6-7
4510	10-11	4515	8-9
4511	12-13	4516	10 & UP
4512	14 & UP		



## NEW! SOFTZONE™

SoftZone products are hand washable, hypoallergenic, and physician recommended. They'll never bottom out and are wear guaranteed for one year.

### WONDERSOLE™

Available in 4 sizes

- Full length insole anatomically designed to provide support to the metatarsals, arch, and heel
- Relieves symptoms related to metatarsalgia, plantar fasciitis, fat pad atrophy, and arthritis
- Soft density colored gel areas provide additional cushioning to bony and sensitive areas

PRODUCT #	MEN'S SIZES	WOMEN'S SIZES
4530	7	5-7
4531	8-9	8-9
4532	10	10-11
4533	11-13	—



### HEEL SPUR WITH REMOVABLE PLUG

- Relief for general or specific heel pain
- Shock-absorbing silicone heel cushion with easily removable plug
- Wear with plug as a comfortable heel pad; wear without the plug for heel relief and protection around the heel spur
- Fits all shoe styles; will not slide in shoe

PRODUCT #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4001	SMALL	—	5-6.5
4002	MEDIUM	5-7.5	7-9.5
4003	LARGE	8-11.5	10-12



### LOW PROFILE HEEL CUSHION

- Low profile fits more shoe styles; will not slide in shoe
- Extra soft silicone heel cushion with shock pattern to cushion and diffuse weight
- Use with the grid down for maximum shock absorption; use with grid up for maximum shear absorption and massaging effect.

PRODUCT #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4004	SMALL	—	5-6.5
4005	MEDIUM	5-7.5	7-9.5
4006	LARGE	8-11.5	10-12



### ¾ LENGTH INSOLE

- Anatomically designed to support longitudinal and transverse arches
- Deeper heel cup that provides superior comfort and control
- Soft, dual-density silicone gel circles provide additional shock absorption and pressure relief to sensitive areas

PRODUCT #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4010	SMALL	—	5-6.5
4011	MEDIUM	5-7.5	7-9.5
4012	LARGE	8-11.5	10-12



### FULL LENGTH INSOLE

- Ideal as a replacement insole in shoes with removable insoles
- Provides relief for heel, arch, and forefoot problems
- Dual formulation absorbs shock under high pressure areas while supporting surrounding tissue
- Multiple vents enhance air circulation

PRODUCT #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4013	SMALL	5-6.5	6.5-8.5
4014	MEDIUM	7-9.5	9-11
4015	LARGE	10-12	11.5-12
4016	X-LARGE	12.5	—



### FULL LENGTH INSOLE WITH TOE CREST

- Thin design ideal for use in shoes with removable insoles or minimal room
- Thin, soft silicone with dual density dots for full-foot comfort and protection
- Broad, soft metatarsal raise for relief of forefoot discomfort
- Designed to reduce pressure while supporting mid- and forefoot areas
- Soft, full toe crest enhances forefoot function

PRODUCT #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4017	SMALL	—	5-6
4018	MEDIUM	5-6	6.5-7.5
4019	LARGE	6.5-7.5	8-9
4033	X-LARGE	8-9	10-10.5
4034	XX-LARGE	9.5-10	11-11.5



### LOW PROFILE FABRIC-COVERED INSOLE

- For general shock absorption and pressure redistribution
- Soft, durable, and comfortable antimicrobial fabric cover allows foot to slide easily into shoe
- May be used under other insoles for shock absorption and protection

PRODUCT #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4030	SMALL	—	4.5-7
4031	MEDIUM	5.5-9.5	7.5-11
4032	LARGE	10-12.5	11.5-12



### GEL METATARSAL PAD WITH TOE SPREADER™

- Anatomically designed to provide comfort and protection of the metatarsal heads
- Helps reduce pain and burning sensations under metatarsal heads
- Ideal for conditions such as Fat Pad Atrophy, Metatarsalgia, Morton's Neuroma, etc.
- Easy to apply; easily fits into dress or casual shoes
- Built-in toe spreader reduces discomfort from Hallux Valgus

PROD. #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4020	SMALL/MEDIUM	5.5-6	4.5-7.5
4021	LARGE/X-LARGE	6.5-11.5	8-11



## 3 Easy Ways to Order

By Phone 800-818-2932

By Fax 908-232-5219

By E-mail [sales@wbcindustries.com](mailto:sales@wbcindustries.com)



### ***WBC Industries: Delivering More™***

Our objective is to help you meet production goals while minimizing your cost, time investment, and on-hand inventory. We're committed to delivering quality orthotic & prosthetic materials at a fair price, with no minimums and no hassles.

At WBC Industries, you're not just another name in our files. We go the extra step to ensure your satisfaction, delivering more than just a product. So whether you need a rush order shipped out right away or you'd just like to chat about your options for custom work, we're there for you.

To find out more about our company and what we can do for you, please visit our web site at [www.wbcindustries.com](http://www.wbcindustries.com) or give us a call at 1-800-818-2932.

**Our full Orthotic & Prosthetic Catalog is also available upon request.**

**Call 1-800-818-2932 to get your free copy!**



625 Central Avenue Westfield NJ 07090

[www.wbcindustries.com](http://www.wbcindustries.com)

1-800-818-2932