

DIGITAL PADS ON A STRIP™

- Elasticized tubing is easy to cut using sharp scissors
- Once cut to size, the tube can be turned inside-out
- End result is virtually the same as Digital Pads™

PROD #	SIZE	USAGE	QUANTITY	SOLD
1105	MINI	5TH TOE	17 PADS/STRIP	2/PKG
1106	S/M	LESSER TOES OR FINGERS	16 PADS/STRIP	2/PKG
1107	L/XL	HALLUX OR THUMB	14 PADS/STRIP	2/PKG

