

FULL LENGTH INSOLE WITH TOE CREST

- Thin design ideal for use in shoes with removable insoles or minimal room
- Thin, soft silicone with dual density dots for full-foot comfort and protection
- Broad, soft metatarsal raise for relief of forefoot discomfort
- Designed to reduce pressure while supporting mid- and forefoot areas
- Soft, full toe crest enhances forefoot function

PRODUCT #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4017	SMALL	—	5-6
4018	MEDIUM	5-6	6.5-7.5
4019	LARGE	6.5-7.5	8-9
4033	X-LARGE	8-9	10-10.5
4034	XX-LARGE	9.5-10	11-11.5

