

LOW PROFILE HEEL CUSHION

- Low profile fits more shoe styles; will not slide in shoe
- Extra soft silicone heel cushion with shock pattern to cushion and diffuse weight
- Use with the grid down for maximum shock absorption; use with grid up for maximum shear absorption and massaging effect.



PRODUCT #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4004	SMALL	—	5-6.5
4005	MEDIUM	5-7.5	7-9.5
4006	LARGE	8-11.5	10-12