

MATERIALS

SUBORTHOLEN®



High strength fleshtone polyethylene that can be cold-formed as well as thermo-molded. Lightweight, yet rigid and strong with minimal creep.

PRODUCT#	SHEET THICKNESS	SIZE
030-0202-00	2.0MM	39" X 39"
030-0203-00	3.0MM	39" X 39"
030-0204-00	4.0MM	39" X 39"
030-0205-00	5.0MM	39" X 39"

NEOPRENE SHEETS



Cost effective generic neoprene available in black in two varieties and three thicknesses. Sheet size is 51" x 83".

PRODUCT#	DESCRIPTION	THICKNESS
15MM/N	NEOPRENE NYLON BOTH SIDES	1.5MM
15MM/UBL	NEOPRENE NYLON UBL/ NYLON	1.5MM
3MM/N	NEOPRENE NYLON BOTH SIDES	3.0MM
3MM/UBL	NEOPRENE NYLON UBL/ NYLON	3.0MM
5MM/N	NEOPRENE NYLON BOTH SIDES	5.0MM
5MM/UBL	NEOPRENE NYLON UBL/ NYLON	5.0MM

FOAM LAMINATE



3/8" foam with 27" of knit loop on both sides. Die cuttable and hook compatible on both sides.

PRODUCT#	COLOR	SIZE
27-FLB	BLACK	YARD

ORTHOPEDIC FELT PADDING



Perfect for casts and pads. 50% Cotton and 50% Wool. Provides abrasion absorbing relief.

PRODUCT#	COLOR	WIDTH	SOLD BY
1WFELT	WHITE	1/8"	YARD
1WFELTW	WHITE	1/4"	YARD

MATERIALS

POLYPROPYLENE PLASTIC SHEETS

PRODUCT#	COLOR	THICKNESS	SHEET SIZE
116PPN2448	NATURAL	1/16"	24" X 48"
116PPN3248	NATURAL	1/16"	32" X 48"
332PPN2448	NATURAL	3/32"	24" X 48"
332PPN3248	NATURAL	3/32"	32" X 48"
18PPN2448	NATURAL	1/8"	24" X 48"
18PPN3248	NATURAL	1/8"	32" X 48"
532PPN2448	NATURAL	5/32"	24" X 48"
532PPN3248	NATURAL	5/32"	32" X 48"
316PPN2448	NATURAL	3/16"	24" X 48"
316PPN3248	NATURAL	3/16"	32" X 48"
14PPN2448	NATURAL	1/4"	24" X 48"
14PPN3248	NATURAL	1/4"	32" X 48"
38PPN2448	NATURAL	3/8"	24" X 48"
38PPN3248	NATURAL	3/8"	32" X 48"
12PPN2448	NATURAL	1/2"	24" X 48"
12PPN3248	NATURAL	1/2"	32" X 48"
116PPW2448	WHITE	1/16"	24" X 48"
116PPW3248	WHITE	1/16"	32" X 48"
332PPW2448	WHITE	3/32"	24" X 48"
332PPW3248	WHITE	3/32"	32" X 48"
18PPW2448	WHITE	1/8"	24" X 48"
18PPW3248	WHITE	1/8"	32" X 48"
532PPW2448	WHITE	5/32"	24" X 48"
532PPW3248	WHITE	5/32"	32" X 48"
316PPW2448	WHITE	3/16"	24" X 48"
316PPW3248	WHITE	3/16"	32" X 48"
14PPW2448	WHITE	1/4"	24" X 48"
14PPW3248	WHITE	1/4"	32" X 48"
18PPRB2448	ROYAL BLUE	1/8"	24" X 48"
18PPRB3248	ROYAL BLUE	1/8"	32" X 48"
532PPRB2448	ROYAL BLUE	5/32"	24" X 48"
532PPRB3248	ROYAL BLUE	5/32"	32" X 48"
316PPRB2448	ROYAL BLUE	3/16"	24" X 48"
316PPRB3248	ROYAL BLUE	3/16"	32" X 48"
18PPL2448	BLACK	1/8"	24" X 48"
18PPL3248	BLACK	1/8"	32" X 48"
532PPL2448	BLACK	5/32"	24" X 48"
532PPL3248	BLACK	5/32"	32" X 48"
316PPL2448	BLACK	3/16"	24" X 48"
316PPL3248	BLACK	3/16"	32" X 48"
14PPL2448	BLACK	1/4"	24" X 48"
14PPL3248	BLACK	1/4"	32" X 48"



Continued →


MATERIALS

POLYPROPYLENE PLASTIC SHEETS, Cont.



PRODUCT#	COLOR	THICKNESS	SHEET SIZE
18PPF2448	FLESHTONE	1/8"	24" X 48"
18PPF3248	FLESHTONE	1/8"	32" X 48"
532PPF2448	FLESHTONE	5/32"	24" X 48"
532PPF3248	FLESHTONE	5/32"	32" X 48"
316PPF2448	FLESHTONE	3/16"	24" X 48"
316PPF3248	FLESHTONE	3/16"	32" X 48"


JMS 500



High molecular weight, high density polyethylene plastic sheet.

PRODUCT#	COLOR	THICKNESS	SHEET SIZE
50001F-1	FLESHTONE	1.0MM	36" X 37.5"
500015F-1	FLESHTONE	1.5MM	36" X 37.5"
50002F-1	FLESHTONE	2.0MM	36" X 37.5"
500025F-1	FLESHTONE	2.5MM	36" X 37.5"
50003F-1	FLESHTONE	3.0MM	36" X 37.5"
500035F-1	FLESHTONE	3.5MM	36" X 37.5"
50004F-1	FLESHTONE	4.0MM	36" X 37.5"
50005F-1	FLESHTONE	5.0MM	36" X 37.5"

CO-POLYMER PLASTIC SHEETS



A blend of polypropylene and 5%-10% low density polyethylene. Sheets are self adhesive when hot and vacuum formable. Co-polymer is less brittle than polypropylene.

PRODUCT#	THICKNESS	SHEET SIZE	COLOR
132361N	3/16"	24" X 48"	NATURAL
132368N	3/16"	32" X 48"	NATURAL
132362N	3/16"	48" X 48"	NATURAL
132363N	3/16"	48" X 96"	NATURAL

MATERIALS

PORON® MEDICAL URETHANE SHEETS

Open-cell urethane is breathable, odorless, washable and non-sensitizing. Will not bottom out over time, and easily bonds to any other orthotic or prosthetic materials.



PRODUCT#	COLOR	THICK	SHEET SIZE	TEXTURE
4708-601	BLUE	1/16"	12" X 42"	SMOOTH-2 SIDES
4708-60118	BLUE	1/16"	18" X 42"	SMOOTH-2 SIDES
4708-60136	BLUE	1/16"	36" X 42"	SMOOTH-2 SIDES
4708-602	BLUE	1/16"	12" X 42"	ABRADED-1 SIDE
4708-60218	BLUE	1/16"	18" X 42"	ABRADED-1 SIDE
4708-60236	BLUE	1/16"	36" X 42"	ABRADED-1 SIDE
4708-603	BLUE	1/16"	12" X 42"	ABRADED-2 SIDES
4708-60318	BLUE	1/16"	18" X 42"	ABRADED-2 SIDES
4708-60336	BLUE	1/16"	36" X 42"	ABRADED-2 SIDES
4708-611	BLUE	1/16"	12" X 42"	ABRADED W/ NAVY CLOTH
4708-61118	BLUE	1/16"	18" X 42"	ABRADED W/ NAVY CLOTH
4708-61136	BLUE	1/16"	36" X 42"	ABRADED W/ NAVY CLOTH
4708-669	BLUE	1/16"	12" X 42"	SCRIM/ ABRADED-1 SIDE
4708-622	TAN	1/16"	12" X 42"	ABRADED-1 SIDE
4708-62218	TAN	1/16"	18" X 42"	ABRADED-1 SIDE
4708-62236	TAN	1/16"	36" X 42"	ABRADED-1 SIDE
4708-623	TAN	1/16"	12" X 42"	ABRADED-2 SIDES
4708-62138	TAN	1/16"	18" X 42"	ABRADED-2 SIDES
4708-62336	TAN	1/16"	36" X 42"	ABRADED-2 SIDES
4708-801	BLUE	1/8"	12" X 42"	SMOOTH-2 SIDES
4708-80118	BLUE	1/8"	18" X 42"	SMOOTH-2 SIDES
4708-80136	BLUE	1/8"	36" X 42"	SMOOTH-2 SIDES
4708-802	BLUE	1/8"	12" X 42"	ABRADED-1 SIDE
4708-80218	BLUE	1/8"	18" X 42"	ABRADED-1 SIDE
4708-80236	BLUE	1/8"	36" X 42"	ABRADED-1 SIDE
4708-803	BLUE	1/8"	12" X 42"	ABRADED-2 SIDES
4708-80318	BLUE	1/8"	18" X 42"	ABRADED-2 SIDES
4708-80336	BLUE	1/8"	36" X 42"	ABRADED-2 SIDES
4708-810	BLUE	1/8"	12" X 42"	ABRADED W/ NAVY CLOTH
4708-81018	BLUE	1/8"	18" X 42"	ABRADED W/ NAVY CLOTH
4708-81036	BLUE	1/8"	36" X 42"	ABRADED W/ NAVY CLOTH
4708-821	TAN	1/8"	12" X 42"	SMOOTH-2 SIDES
4708-8218	TAN	1/8"	18" X 42"	SMOOTH-2 SIDES
4708-82136	TAN	1/8"	36" X 42"	SMOOTH-2 SIDES
4708-822	TAN	1/8"	12" X 42"	ABRADED-1 SIDE
4708-82218	TAN	1/8"	18" X 42"	ABRADED-1 SIDE
4708-82236	TAN	1/8"	36" X 42"	ABRADED-1 SIDE
4708-823	TAN	1/8"	12" X 42"	ABRADED-2 SIDES
4708-82318	TAN	1/8"	18" X 42"	ABRADED-2 SIDES
4708-82336	TAN	1/8"	36" X 42"	ABRADED-2 SIDES

Continued →

MATERIALS

PORON® MEDICAL URETHANE SHEETS, Cont.

PRODUCT#	COLOR	THICK	SHEET SIZE	TEXTURE
4708-869	BLUE	1/8"	12" X 42"	SCRIM/ ABRADED-1 SIDE
4708-86918	BLUE	1/8"	18" X 42"	SCRIM/ ABRADED-1 SIDE
4708-86936	BLUE	1/8"	36" X 42"	SCRIM/ ABRADED-1 SIDE
4708-401	BLUE	1/4"	12" X 42"	SMOOTH-2 SIDES
4708-40118	BLUE	1/4"	18" X 42"	SMOOTH-2 SIDES
4708-40136	BLUE	1/4"	36" X 42"	SMOOTH-2 SIDES
4708-402	BLUE	1/4"	12" X 42"	ABRADED-1 SIDE
4708-40218	BLUE	1/4"	18" X 42"	ABRADED-1 SIDE
4708-40236	BLUE	1/4"	36" X 42"	ABRADED-1 SIDE
4708-403	BLUE	1/4"	12" X 42"	ABRADED-2 SIDES
4708-40318	BLUE	1/4"	18" X 42"	ABRADED-2 SIDES
4708-40336	BLUE	1/4"	36" X 42"	ABRADED-2 SIDES
4708-422	TAN	1/4"	12" X 42"	ABRADED-1 SIDE
4708-42218	TAN	1/4"	18" X 42"	ABRADED-1 SIDE
4708-42236	TAN	1/4"	36" X 42"	ABRADED-1 SIDE
4708-423	TAN	1/4"	12" X 42"	ABRADED-2 SIDES
4708-42318	TAN	1/4"	18" X 42"	ABRADED-2 SIDES
4708-42336	TAN	1/4"	36" X 42"	ABRADED-2 SIDES
4708-301	BLUE	3/8"	12" X 42"	SMOOTH-2 SIDES
4708-30118	BLUE	3/8"	18" X 42"	SMOOTH-2 SIDES
4708-30136	BLUE	3/8"	36" X 42"	SMOOTH-2 SIDES
4708-302	BLUE	3/8"	12" X 42"	ABRADED-1 SIDE
4708-30218	BLUE	3/8"	18" X 42"	ABRADED-1 SIDE
4708-30236	BLUE	3/8"	36" X 42"	ABRADED-1 SIDE
4708-303	BLUE	3/8"	12" X 42"	ABRADED-2 SIDES
4708-30318	BLUE	3/8"	18" X 42"	ABRADED-2 SIDES
4708-30336	BLUE	3/8"	36" X 42"	ABRADED-2 SIDES
4708-321	TAN	3/8"	12" X 42"	SMOOTH-2 SIDES
4708-32118	TAN	3/8"	18" X 42"	SMOOTH-2 SIDES
4708-32136	TAN	3/8"	36" X 42"	SMOOTH-2 SIDES
4708-322	TAN	3/8"	12" X 42"	ABRADED-1 SIDE
4708-32218	TAN	3/8"	18" X 42"	ABRADED-1 SIDE
4708-32236	TAN	3/8"	36" X 42"	ABRADED-1 SIDE
4708-201	BLUE	1/2"	12" X 42"	SMOOTH-2 SIDES
4708-20118	BLUE	1/2"	18" X 42"	SMOOTH-2 SIDES
4708-20136	BLUE	1/2"	36" X 42"	SMOOTH-2 SIDES



MATERIALS

PLASTAZOTE®

Expanded foam by the use of pure nitrogen. Very comfortable and lightweight. Great for orthotics.



PLASTAZOTE® CASE QUANTITIES

- 1" = 3 SHEETS
- 3/4" = 4 SHEETS
- 1/2" = 6 SHEETS
- 3/8" = 8 SHEETS
- 1/4" = 12 SHEETS
- 3/16" = 16 SHEETS
- 1/8" = 24 SHEETS
- 1/16" = 24 SHEETS

Sorry, WBC cannot mix or match thickness or colors to get case prices

PRODUCT#	THICK	SHEET SIZE	COLOR	GRADE & SIZE
705111-1	1/16"	30" X 40"	BLACK	#1 MEDIUM (15 SHORE DUROMETER)
705111-2	1/8"	30" X 40"	BLACK	#1 MEDIUM (15 SHORE DUROMETER)
705111-3	3/16"	30" X 40"	BLACK	#1 MEDIUM (15 SHORE DUROMETER)
705111-4	1/4"	30" X 40"	BLACK	#1 MEDIUM (15 SHORE DUROMETER)
705111-5	3/8"	30" X 40"	BLACK	#1 MEDIUM (15 SHORE DUROMETER)
705111-6	1/2"	30" X 40"	BLACK	#1 MEDIUM (15 SHORE DUROMETER)
705161-2	1/8"	30" X 40"	BLUE	#1 MEDIUM (15 SHORE DUROMETER)
705161-3	3/16"	30" X 40"	BLUE	#1 MEDIUM (15 SHORE DUROMETER)
705161-4	1/4"	30" X 40"	BLUE	#1 MEDIUM (15 SHORE DUROMETER)
705161-6	1/2"	30" X 40"	BLUE	#1 MEDIUM (15 SHORE DUROMETER)
705131-1	1/16"	30" X 40"	WHITE	#1 MEDIUM (15 SHORE DUROMETER)
705131-2	1/8"	30" X 40"	WHITE	#1 MEDIUM (15 SHORE DUROMETER)
705131-3	3/16"	30" X 40"	WHITE	#1 MEDIUM (15 SHORE DUROMETER)
705131-4	1/4"	30" X 40"	WHITE	#1 MEDIUM (15 SHORE DUROMETER)
705131-5	3/8"	30" X 40"	WHITE	#1 MEDIUM (15 SHORE DUROMETER)
705131-6	1/2"	30" X 40"	WHITE	#1 MEDIUM (15 SHORE DUROMETER)
705131-7	3/4"	30" X 40"	WHITE	#1 MEDIUM (15 SHORE DUROMETER)
705131-8	1"	30" X 40"	WHITE	#1 MEDIUM (15 SHORE DUROMETER)
705141-1	1/16"	30" X 40"	PINK	#1 MEDIUM (15 SHORE DUROMETER)
705141-2	1/8"	30" X 40"	PINK	#1 MEDIUM (15 SHORE DUROMETER)
705141-3	3/16"	30" X 40"	PINK	#1 MEDIUM (15 SHORE DUROMETER)
705141-4	1/4"	30" X 40"	PINK	#1 MEDIUM (15 SHORE DUROMETER)
705141-5	3/8"	30" X 40"	PINK	#1 MEDIUM (15 SHORE DUROMETER)
705141-6	1/2"	30" X 40"	PINK	#1 MEDIUM (15 SHORE DUROMETER)
705141-7	3/4"	30" X 40"	PINK	#1 MEDIUM (15 SHORE DUROMETER)
705141-8	1"	30" X 40"	PINK	#1 MEDIUM (15 SHORE DUROMETER)
705331-1	1/16"	34" X 34"	WHITE	#2 FIRM (SHORE A DUROMETER)
705331-2	1/8"	34" X 34"	WHITE	#2 FIRM (SHORE A DUROMETER)
705331-3	3/16"	34" X 34"	WHITE	#2 FIRM (SHORE A DUROMETER)
705331-4	1/4"	34" X 34"	WHITE	#2 FIRM (SHORE A DUROMETER)
705331-5	3/8"	34" X 34"	WHITE	#2 FIRM (SHORE A DUROMETER)
705331-6	1/2"	34" X 34"	WHITE	#2 FIRM (SHORE A DUROMETER)
705331-7	3/4"	34" X 34"	WHITE	#2 FIRM (SHORE A DUROMETER)
705411-1	1/16"	25" X 37"	BLACK	#3 RIGID (60 SHORE A DUROMETER)
705411-2	1/8"	25" X 37"	BLACK	#3 RIGID (60 SHORE A DUROMETER)
705411-3	3/16"	25" X 37"	BLACK	#3 RIGID (60 SHORE A DUROMETER)
705411-4	1/4"	25" X 37"	BLACK	#3 RIGID (60 SHORE A DUROMETER)
705411-5	3/8"	25" X 37"	BLACK	#3 RIGID (60 SHORE A DUROMETER)
705411-6	1/2"	25" X 37"	BLACK	#3 RIGID (60 SHORE A DUROMETER)

MATERIALS

PELITE®

Expanded lightweight closed-cell polyethylene foam.

PROD#	THICKNESS	SIZE	COLOR	DENSITY
1103-01	1/8" (3.0MM)	39" X 39"	WHITE	FIRM DENSITY, NON-PERFORATED
1105-00	3/16" (5.0MM)	39" X 39"	WHITE	FIRM DENSITY, NON-PERFORATED
1108-01	5/16" (8.0MM)	39" X 39"	WHITE	FIRM DENSITY, NON-PERFORATED
1110-01	3/8" (10.0MM)	39" X 39"	WHITE	FIRM DENSITY, NON-PERFORATED
1113-01	1/2" (13.0MM)	39" X 39"	WHITE	FIRM DENSITY, NON-PERFORATED
1306-06	1/8" (3.0MM)	39" X 39"	WHITE	MEDIUM DENSITY, NON-PERFORATED
1313-16	1/8" (3.0MM)	39" X 39"	WHITE	MEDIUM DENSITY, NON-PERFORATED

P-CELL™ SHEETS

New heat moldable, grindable closed-cell structure. Superior in cushioning, shock absorption, tensile strength, moisture absorption, and abrasiveness.

PRODUCT#	THICKNESS	SHEET SIZE	COLOR
SDUR-00159	1/16"	44" X 30"	PINK
SDUR-00161	1/8"	44" X 30"	PINK
SDUR-00163	3/16"	44" X 30"	PINK
SDUR-00165	1/4"	44" X 30"	PINK
SDUR-00167	1/2"	44" X 30"	PINK
SDUR-00169	3/4"	44" X 30"	PINK
SDUR-00171	1"	44" X 30"	PINK


NEOLON (AKA SPENCO)

Our new and improved firmer neoprene - for more active patients. Colored nylon top cover with bottom side abraded.

PRODUCT#	THICKNESS	COVER COLOR	SIZE
416	1/16"	ROYAL BLUE	41" X 51"
426	1/16"	ROYAL BLUE	82" X 51"
418	1/8"	ROYAL BLUE	41" X 51"
428	1/8"	ROYAL BLUE	82" X 51"

MATERIALS

NEOLON, Cont.



PRODUCT#	THICKNESS	COVER COLOR	SIZE
516	1/16"	BLACK	41" X 51"
526	1/16"	BLACK	82" X 51"
518	1/8"	BLACK	41" X 51"
528	1/8"	BLACK	82" X 51"
616	1/16"	WINE	41" X 51"
626	1/16"	WINE	82" X 51"
618	1/8"	WINE	41" X 51"
628	1/8"	WINE	82" X 51"
718	1/8"	GREEN	41" X 51"
728	1/8"	GREEN	82" X 51"


WBC EVA FOAM SHEETS



Used for foot orthoses and prosthetic devices, orthotic posting and heel lifts. Can be used as a base, middle layer or top cover for many foot orthoses. Highly resilient and washable.

PRODUCT#	COLOR	THICKNESS	DUROMETER	SHEET SIZE
169-10-02	PURPLE SWIRL	2.0MM	35	25" X 43"
169-10-03	PURPLE SWIRL	3.0MM	35	25" X 43"
169-20-02	RED	2.0MM	35	25" X 43"
169-20-03	RED	3.0MM	35	25" X 43"
169-30-02	WATERCOLOR	2.0MM	35	25" X 43"
169-30-03	WATERCOLOR	3.0MM	35	25" X 43"

COMFORT CELL





Combination blended material of EVA and a touch of neoprene producing a better rebound, longer lasting material. It molds and grinds easily, and produces lower dust than straight EVA. Possesses superb cushioning and shape retention. Available in a variety of colors, durometers, thicknesses, and combinations.

PRODUCT#	DIMENSIONS	THICKNESS	DUROMETER	COLOR
30001-3	40" X 36"	1/16", 1/8", 1/4"	30	WHITE
30011-13	40" X 36"	1/16", 1/8", 1/4"	30	POWDER BLUE
30021-23	40" X 36"	1/16", 1/8", 1/4"	30	PINK


Continued →

MATERIALS



COMFORT CELL, Cont.

PROD#	DIMENSIONS	THICKNESS	DURO.	COLOR
40001-3	40" X 36"	1/8", 1/4", 3/8"	40	WHITE
40011-15	40" X 36"	1/16", 1/8", 3/16", 1/4", 3/8", 1/2"	40	ROYAL BLUE
40021-23	40" X 36"	1/16", 1/8", 1/4"	40	RED
40031-34	40" X 36"	1/16", 1/8", 1/4", 3/8", 1/2"	40	DENIM BLUE
40041-43	40" X 36"	1/16", 1/8", 1/4"	40	TEAL
35001-3	40" X 36"	1/16", 1/8", 1/4"	35	GOLD, WHITE, BLACK
35011-13	40" X 36"	1/16", 1/8", 1/4"	35	RED, WHITE, BLACK
35021-23	40" X 36"	1/16", 1/8", 1/4"	35	TEAL, WHITE, BLACK
50001-5	38" X 32"	1/16", 1/8", 1/4", 3/8", 1/2"	50	BLACK
60001-5	36" X 32"	1/16", 1/8", 1/4", 3/8", 1/2"	60	WHITE
60011-15	36" X 32"	1/8", 3/16", 1/4", 3/8", 1/2"	60	TEAL
80001-8	27" X 31"	1/16", 1/8", 3/16", 1/4", 3/8", 1/2", 5/8", 3/4"	80	BLACK
80011-14	27" X 31"	1/8", 1/4", 1/2", 3/4"	80	WHITE
80021-25	27" X 31"	1/8", 3/16", 1/4", 1/2", 3/4"	80	BROWN
80031-34	27" X 31"	1/8", 1/4", 1/2", 3/4"	80	BEIGE



COMFORT CELL DUAL LAMINATIONS

Produced with softer durometer top and firmer durometer bottom. Perfect for orthotic foot beds and general padding. Has a heating temperature of 250° - 275° F and is easily grindable. Comes in pink and royal blue, or white and royal blue.

PROD#	DIMENSIONS	THICKNESS	TOP LAYER	BOTTOM LAYER
22003	38" X 35"	1/4"	1/8" PINK 30 DURO	1/8" ROYAL BLUE 40 DURO
22204	38" X 36"	1/4"	1/8" WHITE 30 DURO	1/8" ROYAL BLUE 40 DURO

COMFORT CELL TRI-LAMINATIONS

Available with Urethane, Ameriprene, or Amerithane center layers. Urethane provides durability and longevity. Ameriprene and Amerithane provides cushioning and moldability, while being a cost effective alternative to urethane center layers.



PRODUCT#	DIMENSIONS	THICKNESS	CENTER MATERIAL
33001-2,4,5-10	38" X 29"	1/4" W/ ASSORTED TOP & BOTTOM	URETHANE
33011-12-13	38" X 29"	1/4" W/ ASSORTED TOP & BOTTOM	AMERIPRENE


Continued →

MATERIALS

COMFORT CELL TRI-LAMINATIONS, Cont.




PRODUCT#	DIMENSIONS	THICKNESS	CENTER MATERIAL
33014-15-16	38" X 32"	1/4" W/ ASSORTED TOP & BOTTOM	AMERITHANE
33017	38" X 29"	5/16" W/ PINK TOP & BLUE BOTTOM	URETHANE
33018	38" X 29"	5/16" W/ PINK TOP & BLUE BOTTOM	AMERIPRENE
33019	38" X 29"	5/16" W/ PINK TOP & BLUE BOTTOM	AMERITHANE



SOLING & TOP COVER MATERIAL

Soling is available in four standard colors. Heel & wedging material is available in matching colors. Soling can be cut in any direction to minimize waste and is produced in a tight waffle pattern. Top cover material helps disperse moisture and provides a four way stretch along with excellent cushioning.

PRODUCT#	DIMENSIONS	THICKNESS	DUROMETER	COLOR
99001	22" X 36"	5/16"	60	BLACK
99002	22" X 36"	5/16"	60	BROWN
99003	22" X 36"	5/16"	60	BEIGE
99004	22" X 36"	5/16"	60	WHITE
10001	36" X 26"	1/16"	AMERIPRENE W/ TOP COVER	BLACK
10002	36" X 26"	1/8"	AMERIPRENE W/ TOP COVER	BLACK



WBC LAMINATED SHEETS


Easy to grind, cut, and/ or shape to any size using standard equipment. Adheres easily to other materials with contact cement. EVA (aka Thermosky®) PPT (aka Poron®).

PROD#	DESCRIPTION	THICK	SHEET SIZE
551	1/8" BLUE EVA & 1/8" PINK PLASTAZOTE	1/4"	40" X 60"
5512	1/8" BLUE EVA & 1/8" PINK PLASTAZOTE	1/4"	30" X 40"
552	1/4" BLUE EVA & 1/8" PINK PLASTAZOTE	3/8"	40" X 60"
5522	1/4" BLUE EVA & 1/8" PINK PLASTAZOTE	3/8"	30" X 40"
533	1/4" BLUE EVA & 1/4" PINK PLASTAZOTE	1/2"	40" X 60"
5532	1/4" BLUE EVA & 1/4" PINK PLASTAZOTE	1/2"	30" X 40"
21110	1/8" BLUE EVA, 1/16" BLUE PORON, 1/8" BLACK EVA	5/16"	30" X 40"
21130	1/16" BLUE EVA, 1/16" BLUE PORON, 1/8" BLUE EVA	1/4"	30" X 40"
21131	1/16" WHITE EVA, 1/16" BLUE PORON, 1/8" PINK EVA	1/4"	30" X 40"

Continued →

MATERIALS


WBC LAMINATED SHEETS, Cont.



PROD#	DESCRIPTION	THICK	SHEET SIZE
21132	1/16" PINK EVA, 1/16" BLUE PORON, 1/8" PINK EVA	1/4"	30" X 40"
21133	1/8" BLUE EVA, 1/8" BLUE FORM EVA, 1/8" BLUE EVA	3/8"	31" X 28"
21135	1/4" PINK PLASTAZOTE, 1/16" BLUE PORON, 1/8" BLUE EVA	7/16"	30" X 40"
21137	1/8" PINK PLASTAZOTE, 1/16" BLUE PORON, 1/8" BLUE EVA	5/16"	30" X 40"

DUAL LAMINATES

P-CELL + PORON



PRODUCT#	DESCRIPTION	THICK	SHEET SIZE
PCPN-00014	1/8" P-CELL + 1/16" BLUE PORON	3/16"	37" X 14"
PCPN-00006	1/8" P-CELL + 1/16" BLUE PORON	3/16"	37" X 28"
PCPN-00015	1/8" P-CELL + 1/8" BLUE PORON	1/4"	37" X 14"
PCPN-00007	1/8" P-CELL + 1/8" BLUE PORON	1/4"	37" X 28"
PCPN-00013	1/4" P-CELL + 1/8" BLUE PORON	3/8"	37" X 14"
PCPN-00005	1/4" P-CELL + 1/8" BLUE PORON	3/8"	37" X 28"


PLASTAZOTE + PORON



PRODUCT#	DESCRIPTION	THICK	SIZE
PZPN-00014	1/8" #1 PINK PLASTAZOTE + 1/16" BLUE PORON	3/16"	37" X 14"
PZPN-00006	1/8" #1 PINK PLASTAZOTE + 1/16" BLUE PORON	3/16"	37" X 28"
PZPN-00015	1/8" #1 PINK PLASTAZOTE + 1/8" BLUE PORON	1/4"	37" X 14"
PZPN-00007	1/8" #1 PINK PLASTAZOTE + 1/8" BLUE PORON	1/4"	37" X 28"
PZPN-00013	1/4" #1 PINK PLASTAZOTE + 1/8" BLUE PORON	3/8"	37" X 14"
PZPN-00005	1/4" #1 PINK PLASTAZOTE + 1/8" BLUE PORON	3/8"	37" X 28"

MATERIALS

TRIPLE LAMINATES



P-CELL + PORON + MICROCEL PUFF

PRODUCT#	DESCRIPTION	THICK	SHEET SIZE
MDUR-00083	1/8" P-CELL + 1/16" BLUE PORON + 1/8" BLUE PUFF	5/16"	37" X 30"
MDUR-00084	1/8" P-CELL + 1/16" BLUE PORON + 1/4" BLUE PUFF	7/16"	37" X 30"
MDUR-00085	1/4" P-CELL + 1/16" BLUE PORON + 1/8" BLUE PUFF	7/16"	37" X 30"
MDUR-00086	1/8" P-CELL + 1/8" BLUE PORON + 1/8" BLUE PUFF	3/8"	37" X 30"

PLASTAZOTE + PORON + MICROCEL PUFF



PRODUCT#	DESCRIPTION	THICK	SHEET SIZE
MDUR-00001	1/8" #1 PINK PLASTAZOTE + 1/16" BLUE PORON + 1/8" BLUE PUFF	5/16"	37" X 30"
MDUR-00002	1/8" #1 PINK PLASTAZOTE + 1/16" BLUE PORON + 1/4" BLUE PUFF	7/16"	37" X 30"
MDUR-00003	1/4" #1 PINK PLASTAZOTE + 1/16" BLUE PORON + 1/8" BLUE PUFF	7/16"	37" X 30"
MDUR-00058	1/8" #1 PINK PLASTAZOTE + 1/8" BLUE PORON + 1/8" BLUE PUFF	3/8"	37" X 30"

TRILAM



PRODUCT#	DESCRIPTION	THICK	SHEET SIZE
TRLM-00001	1/16" PINK PUFF + 1/16" BLUE PORON + 1/8" PINK PUFF	1/4"	37" X 30"
TRLM-00007	1/16" BBA PUFF + 1/16" BLUE PORON + 1/8" PINK PUFF	1/4"	37" X 30"
TRLM-00003	1/16" BLUE PUFF + 1/16" BLUE PORON + 1/8" BLUE PUFF	1/4"	37" X 30"
TRLM-00005	1/16" WHITE PUFF + 1/16" BLUE PORON + 1/8" WHITE PUFF	1/4"	37" X 30"

MATERIALS

VOLARA 4E

A soft density, closed-cell, cross linked polyethylene foam padding material. It is easy to mold and sand, auto-adhesive and great for cushioning orthoses.

PRODUCT#	SIZE	COLOR
4 E125R	1/8" X 30" X 50 YARDS	WHITE
4 E187R	3/16" X 30" X 50 YARDS	WHITE
4 E250R	1/4" X 30" X 50 YARDS	WHITE

VOLARA 6A

Medium density. Thermoplastic heat moldable foam.

PRODUCT#	SIZE	COLOR
6 A125R	1/8" X 30" X 50 YARDS	WHITE

SURLYN®

Transparent, long lasting with durable wear. Outstanding thermo forming characteristics.

PRODUCT#	THICKNESS	SHEET SIZE
SL14007	1/8"	24" X 48"
SL14008	1/8"	32" X 48"
SL14004	1/8"	36" X 48"
SL14017	3/16"	24" X 48"
SL14018	3/16"	32" X 48"
SL14014	3/16"	36" X 48"
SL14027	1/4"	24" X 48"
SL14028	1/4"	32" X 48"
SL14024	1/4"	36" X 48"
SL13996	3/8"	16" X 16"
SL13995	1/2"	16" X 16"

MATERIALS

VIVAK®

AKA Durr-Plex, Diaclean, PETG, Thermocheck


Transparent color great for check sockets (rigid, self adhesive when hot, vacuum formable).



PRODUCT#	THICKNESS	SHEET SIZE
129181	3.0MM	12" X 12"
129182	3.0MM	16" X 16"
129183	3.0MM	24" X 48"
129361	4.5MM	12" X 12"
129362	4.5MM	16" X 16"
129363	4.5MM	24" X 48"
129251	6.0MM	12" X 12"
129252	6.0MM	16" X 16"
129253	6.0MM	24" X 48"
129371	9.0MM	12" X 12"
129372	9.0MM	16" X 16"
129373	9.0MM	24" X 48"

LOW DENSITY POLYETHYLENE (LDPE)

Soft and flexible material, ideal for body jackets. Excellent flexibility with good formability.




PRODUCT#	THICKNESS	SHEET SIZE	COLOR
145161N	1/16"	24" X 48"	NATURAL
145168N	1/16"	32" X 48"	NATURAL
145162N	1/16"	48" X 48"	NATURAL
145163N	1/16"	48" X 96"	NATURAL
145181N	1/8"	24" X 48"	NATURAL
145188N	1/8"	32" X 48"	NATURAL
145182N	1/8"	48" X 48"	NATURAL
145183N	1/8"	48" X 96"	NATURAL
145531N	5/32"	24" X 48"	NATURAL
145538N	5/32"	32" X 48"	NATURAL
145532N	5/32"	48" X 48"	NATURAL
145533N	5/32"	48" X 96"	NATURAL

Continued →


MATERIALS

LOW DENSITY POLYETHYLENE (LDPE), Cont.



PRODUCT#	THICKNESS	SHEET SIZE	COLOR
145361N	3/16"	24" X 48"	NATURAL
145368N	3/16"	32" X 48"	NATURAL
145362N	3/16"	48" X 48"	NATURAL
145363N	3/16"	48" X 96"	NATURAL
145251N	1/4"	24" X 48"	NATURAL
145258N	1/4"	32" X 48"	NATURAL
145252N	1/4"	48" X 48"	NATURAL
145253N	1/4"	48" X 96"	NATURAL
145121N	1/2"	24" X 48"	NATURAL
145122N	1/2"	48" X 48"	NATURAL
145123N	1/2"	48" X 96"	NATURAL

HIGH DENSITY POLYETHYLENE (HDPE)



Very durable product and very good thermo forming, but not recommended for dynamic orthosis applications requiring constant flexibility.

PRODUCT#	THICKNESS	SHEET SIZE	COLOR
140161N	1/16"	24" X 48"	NATURAL
140162N	1/16"	48" X 48"	NATURAL
140163N	1/16"	48" X 96"	NATURAL
140331N	3/32"	24" X 48"	NATURAL
140332N	3/32"	48" X 48"	NATURAL
140333N	3/32"	48" X 96"	NATURAL
140181N	1/8"	24" X 48"	NATURAL
140182N	1/8"	48" X 48"	NATURAL
140183N	1/8"	48" X 96"	NATURAL
140531N	5/32"	24" X 48"	NATURAL
140532N	5/32"	48" X 48"	NATURAL
140533N	5/32"	48" X 96"	NATURAL
140361N	3/16"	24" X 48"	NATURAL
140362N	3/16"	48" X 48"	NATURAL
140363N	3/16"	48" X 96"	NATURAL

MATERIALS

ORTHOTIC CORK MATERIALS

8012 CUSHION CORK



A cork and rubber mix that assures flexibility and durability. Easy to grind. Not heat moldable.

PRODUCT#	DESCRIPTION	THICKNESS	SHEET SIZE
8012-1	CUSHION CORK	1/16"	18" X 36"
8012-2	CUSHION CORK	1/8"	18" X 36"
8012-3	CUSHION CORK	3/16"	18" X 36"
8012-4	CUSHION CORK	1/4"	18" X 36"
8012-5	CUSHION CORK	3/8"	18" X 36"
8012-6	CUSHION CORK	1/2"	18" X 36"


8013 THERMOKORK



A thermoplastic orthotic cork that is heat formable to 170°F. A technician can form the cork directly to the positive mold. Thermokork is strong, flexible, and will not "pack" or "bottom." Since it is light weight, it is tolerated by most patients. Cover with Plastazote®, S.T.S., or leather to customize your orthoses to the appropriate medical condition. Easy to grind and trim.

PRODUCT#	DESCRIPTION	THICKNESS	SHEET SIZE
8013-2	THERMOKORK	1/8"	36" X 44"
8013-3	THERMOKORK	3/16"	36" X 44"
8013-4	THERMOKORK	1/4"	36" X 44"
8013-5	THERMOKORK	3/8"	36" X 44"

8014 RUBBER CORK



A simulated cork which contains rubber. Used for rear foot posting and as a filler under the arch. Not heat moldable.

PRODUCT#	DESCRIPTION	THICKNESS	SHEET SIZE
8014-1	RUBBER CORK	1/16"	18" X 36"
8014-2	RUBBER CORK	1/8"	18" X 36"
8014-3	RUBBER CORK	3/16"	18" X 36"
8014-4	RUBBER CORK	1/4"	18" X 36"
8014-5	RUBBER CORK	3/8"	18" X 36"

MATERIALS

BIRKO-CORK



Made from natural cork granules mixed with nylon. Birko-Cork is porous, lightweight, flexible, and thermopliable. Material becomes pliable, but not sticky at 174°F. At 212°F it is slightly sticky and good for molding. At 250°F it is sticky, excellent for molding, and adheres well to other materials. It is also economical in use, since leftovers from job can be kneaded together for future use.

PROD#	DESCRIPTION	THICKNESS	SHEET SIZE	DENSITY
13105	BIRKO-CORK	3/16" (5.0MM)	12" X 18"	SOFT
13010	BIRKO-CORK	3/8" (10.0MM)	12" X 18"	SOFT
13305	BIRKO-CORK	3/16" (5.0MM)	12" X 18"	HARD
13310	BIRKO-CORK	3/8" (10.0MM)	12" X 18"	HARD


WBC BIO-KORK



A traditional blend of EVA & cork materials. WBC cork is 60-65 durometer, vacuum formable and easy to grind. Heat formable at 300°F.

PROD#	DESCRIPTION	THICKNESS	SHEET SIZE
33003	WBC BIO-KORK	1/8" (3.0MM)	32" X 40"
33004	WBC BIO-KORK	5/32" (4.0MM)	32" X 40"
33006	WBC BIO-KORK	1/4" (8.0MM)	32" X 40"
33008	WBC BIO-KORK	5/16" (8.0MM)	32" X 40"
33010	WBC BIO-KORK	3/8" (10.0MM)	32" X 40"
33012	WBC BIO-KORK	1/2" (12.0MM)	32" X 40"

AMERICORK



A tan colored imitation cork material giving moldability and durability along with light weight formulation. Ideal to replace cork for most applications. Use for build ups, reinforcements, orthotic base layers, and postings. Grindable and has a heating temperature of 275° - 300° F.

PROD#	DIMENSIONS	THICKNESS	DUROMETER
65001-4	35" X 32"	1/8", 1/4", 3/8", 1/2"	60 - 65