

SILIPOS® BRAND PRODUCTS

DIGITAL PAD™



Protects toes and fingers from friction and toes from footwear pressure. Use where corns, calluses, or hammer toes cause discomfort.

PROD #	SIZE	USAGE	SOLD
10055	MINI	5TH TOE	12/ PKG
10065	S/M	LESSER TOES OR FINGERS	12/ PKG
10075	L/ XL	HALLUX OR THUMB	12/ PKG



DIGITAL PADS ON A STRIP™

Pads are on a continuous strip. Elasticized tubing is easy to cut using sharp scissors. Once cut to size, the tube can be turned inside-out. The end result is virtually the same as the Digital Pad (above). Savings on labor cost make the strip very cost effective.

PROD #	SIZE	USAGE	QTY	SOLD
1105	MINI	5TH TOE	17 PADS/ STRIP	2/ PKG
1106	S/M	LESSER TOES OR FINGERS	16 PADS/ STRIP	2/ PKG
1107	L/ XL	HALLUX OR THUMB	14 PADS/ STRIP	2/ PKG



GEL TUBING™

Fully lined with Silipos® Polymer Gel. Features all-round digital protection. Uses less room than conventional foam tubing. Gel tubing is durable and effective. The gel will not bottom out as foam materials do. Can help reduce scar tissue and can be cut to size required for toes or fingers.

PROD #	STYLE	SIZE	SOLD
10865	RIBBED KNIT FABRIC	NARROW, 3/4" X 6"	2/PKG
10875	RIBBED KNIT FABRIC	WIDE, 1" X 6"	2/PKG
10855	MESH	NARROW, 3/4" X 6"	3/PKG
10845	MESH	WIDE, 3" X 10"	1/PKG

SILIPOS® BRAND PRODUCTS



GEL TOE CREST

100% medical grade silicone absorbs shock & vibration and provides even pressure distribution. Helps reduce pressure and pain caused by hammer, claw, and mallet toes. Will not bottom out as foam materials do.

PRODUCT #	SIZE	STYLE	SOLD
10415	SMALL	LEFTS	3 IN A POLY BAG
10425	SMALL	RIGHTS	3 IN A POLY BAG
10435	LARGE	LEFTS	3 IN A POLY BAG
10445	LARGE	RIGHTS	3 IN A POLY BAG
10515	SMALL	LEFT	1 IN A BOX
10525	SMALL	RIGHT	1 IN A BOX
10535	LARGE	LEFT	1 IN A BOX
10545	LARGE	RIGHT	1 IN A BOX



DIGITAL CAPS™

Fully lined with Silipos® Polymer Gel, the Digital Cap is ideal for all-round digital protection. Widely used in sports for pressure and friction protection. Can be used to protect over/ under lapping toes. Available in two materials - stretchable fabric material offering a mild compression depending on sizes chosen, and mesh material that is more elasticized. They are long enough for fingers and fabric will not fray when cut. Average cap length is 2 1/2".

PRODUCT #	STYLE	SIZE	SOLD
10195	RIBBED KNIT FABRIC	MINI	6/ PKG
10205	RIBBED KNIT FABRIC	S/M	6/ PKG
10206	RIBBED KNIT FABRIC	L/XL	6/ PKG
10235	MESH	ONE SIZE	8/ PKG

SILIPOS® BRAND PRODUCTS

TOE SEPARATORS™



Anatomically molded to fit comfortably between toes reducing shearing and inter-digital pressure. Can be trimmed if required and is a cost-effective alternative to silicone compounds.

PRODUCT #	SIZE	SOLD
11605	SMALL	15/ BAG
11615	LARGE	15/ BAG

TOE SPREADERS™



Thicker than the separators (above) due to the polymer's soft conforming properties, these spreaders are extremely comfortable used as an inter-digital wedge for bunion deformity, overlapping toes or toe drift. They can be trimmed with scissors.

PRODUCT #	SIZE	SOLD
11505	SMALL	4/ BAG
11515	MEDIUM	4/ BAG
11525	LARGE	4/ BAG

DIGISTRIP™



Half coated with Silipos® Polymer Gel. Very cost effective and approximately 24" (60 cm) long. Same applications as the digital pad but in a strip, allowing you to choose the length needed. Working on 3 cm per toe, it will yield 20 toe protectors per box.

PRODUCT #	SIZE	SOLD
1100	NARROW	2 STRIPS/ PKG
1101	MEDIUM	2 STRIPS/ PKG
1102	WIDE	2 STRIPS/ PKG

SILIPOS® BRAND PRODUCTS

ALL GEL DIGITAL CAP



Ideal for easing pressure and friction caused by corns, blisters, calluses, hammertoes, and in-grown nails. Fully lined with Silipos® gel to moisturize, condition, and soften the skin. Virtually transparent and perfect for open toe shoes/ sandals. Average cap length is 2.5"; and they are washable and reusable.

PRODUCT #	SIZE	SOLD
10495	SMALL/ MEDIUM	2/ PACK
10505	LARGE/ X-LARGE	2/ PACK



ARTHRITIC/ DIABETIC GEL SOCK™

Single layer sock with 1/16" (2mm) of gel throughout the sole. Helps prevent foot ulcerations. Features a seamless toe and non restrictive cuff. May be used with or without diabetic shoes. Provides optimum protection for the neuropathic foot. Helps reduce friction, abrasion, shear forces & callus formation. Serves as a soft tissue supplement for fat pad atrophy. One pair per package.

PRODUCT #	U.S.	OTHER			
WHITE	BLACK	SOCK SIZE	SOCK SIZE	CM	SHOE SIZE
1701	1706	7-9	37-40	25-27	FEMALE SIZE 6 OR UNDER
1702	1707	9-11	41-43	27-29	FEMALE OVER SIZE 6/ MALE UNDER SIZE 9
1703	1708	10-13	44-46	28-31	MALE SIZE 9 AND OVER

SILIPOS® BRAND PRODUCTS



ACHILLES HEEL PAD™

A soft conforming gel pad protects the posterior of the heel from pressure and friction. Provides gentle compression and dissipation of pressure around an inflamed area. Effective way to protect back of heel from pressure and friction caused by footwear pressure from sports boots: ski, in-line, hockey, ice skating. Also, reduces pain from footwear pressure particularly on any irregular heel bumps.

(Gel Thickness - 1/8")

PRODUCT #	SIZE	SOLD
10385	S/M	1/ PKG
10395	L/XL	1/PKG



MALLEOLAR SLEEVE™

Two large discs of polymer gel molded on both medial and lateral malleolar, reducing pressure and friction to the area. Can be used pre or post injury. Provides gentle compression to ankle bones and an effective way to protect ankles from sports boots: ski, in-line, hockey, and ice skating. Reduces abrasion and shear forces.

PRODUCT #	SIZE	SOLD
15005	ONE SIZE	1/ PKG



HEEL/ ELBOW SLIP OVER™

FOR HEEL: Helps the healing of dry, cracked heels or elbows. Used to help prevent heel ulceration. Gel pad can be placed on dorsal side of foot to absorb shoe lace bite, or can be positioned anywhere around the foot/ heel. Features protective liner for orthopedic braces or an AFO.

FOR ELBOW: Protective padding for ulnar nerve.

PRODUCT #	SIZE	SOLD
15225	S/M	1/ PKG
15235	L/XL	1/ PKG

SILIPOS® BRAND PRODUCTS

SOFT SKIN HEEL SLEEVE™



Gel pad anatomically molded to the shape of the heel. Will moisturize and soften dry, hard skin, helping to reduce heel fissures. Occlusion has been shown to be beneficial for psoriasis and eczema. Can be used in conjunction with skin moisturizing cream.

PRODUCT #	SIZE	SOLD
15265	ONE SIZE	1/ PKG

UNIVERSAL GEL STRAP™



Surgical grade elastic sleeve with a gel pad that covers all five metatarsal heads. Helps reduce pain and discomfort. Serves as a soft tissue supplement to fat pad atrophy. Pad is available uncovered for maximum benefit of the mineral oil based gel helping to soften and reduce callused skin. Pad is also available covered to provide an extended life when callusing is not a problem. By cutting the stitch, it can be also used on the instep, elbow, wrist, or hand.

PRODUCT #	SIZE	SOLD
10155	S/M (GEL UNCOVERED)	1/ PKG
10165	L/XL (GEL UNCOVERED)	1/ PKG
10175	S/M (GEL COVERED)	1/ PKG
10185	L/XL (GEL COVERED)	1/ PKG

BALL OF FOOT GEL CUSHION™



Excellent for relieving the burning sensation under the metatarsal heads of the foot. Protects against shear forces, abrasion, friction, and pressure. May be used as a soft tissue supplement for metatarsalgia. Ideal for calluses and small neuromas. Gel thickness is 2mm.

PRODUCT #	SOLD
10025	2 PCS./ PKG

SILIPOS® BRAND PRODUCTS

GEL METATARSAL™



Designed to protect the metatarsals. Reduces friction and pressure. Flexible elastic loop can be worn over either the second or third toe. Can be worn with casual/ dress shoes.

PRODUCT #	SIZE	SOLD
10125	SMALL	1/ PKG
10135	LARGE	1/ PKG

GEL METATARSAL PAD™



100% medical grade silicone absorbs shock, vibration, and provides even pressure distribution. Anatomically designed for comfort and support. Helps reduce the burning sensation under your metatarsal heads. Ideal for conditions such as fat pad atrophy, metatarsalgia, morton's neuroma, etc.

PRODUCT #	SIZE	SOLD
10465	ONE SIZE	1 PAIR/ PKG

GEL FOOT COVER™



Made from a thin layer of gel covered on both sides with fabric. Offers complete toe protection and extends to under the metatarsal heads. Worn comfortably in everyday footwear. Ideal if activity places stress on toes. Can be worn under hosiery.

PROD #	SIZE	SOLD
10275	SMALL - 3.75" X 4.9" (9.5 cm X 12.5 cm)	1 PAIR/ PKG
10285	MEDIUM - 4.0" X 5.3" (10 cm X 13.5 cm)	1 PAIR/ PKG
10295	LARGE - 4.3" X 5.9" (11 cm X 15 cm)	1 PAIR/ PKG

SILIPOS® BRAND PRODUCTS



BUNIONCARE GEL SLEEVE™

Comfortable sleeve cushions, pads, and protects. Anatomically designed gel pad conforms to the shape of the bunion. Alleviates friction and pressure. Helps reduce scarring after surgery.

PRODUCT #	SIZE	SOLD
10325	S/M	1 PAIR/ PKG
10335	L/XL	1 PAIR/ PKG



PROTECTIVE BUNION SHIELD™

Cushions and conforms to painful bone protrusion. Protects bunion joint and the hallux. Fully lined with Silipos® Polymer Gel. Fabric can be cut between 1st & 2nd digit to allow correct positioning. Elastic strap holds appliance in place.

PRODUCT #	SIZE	SOLD
10105	S/M	4/ PKG
10115	L/XL	4/ PKG



GEL BUNION™

Protects bunion joint using a concave gel pad. Soft foam shell for additional pressure absorption. Toe loop can be cut at the spoon end to give comfortable fit for each patient. Gel pad comes separate from shell, allowing exact positioning for each patient. Can be worn with dress shoes.

PRODUCT #	SIZE	SOLD
10085	SMALL	5/ PKG
10095	LARGE	5/ PKG

SILIPOS® BRAND PRODUCTS



BOOT BUMPER™

Conforming polymer gel pads that protect the front and back of the ankle from pressure, pinching, and friction. Helps protect the foot and ankle from abrasions caused by athletic footwear, without the added bulk. Effective way to protect ankles from sports boots: ski, in-line, hockey, and ice skating.


PRODUCT #	SIZE	SOLD
1770	S/M	1/ PKG
1772	L/XL	1/ PKG



WONDERSPUR™

Available in 4 sizes, effective in treating heel spurs; plantar fasciitis; knee, back & ankle pain; and fat pad atrophy.

PROD #	SIZE	U.S. MEN'S SIZES	U.S. WOMEN'S SIZES
4500	SMALL	2-4	4-5
4501	MEDIUM	5-7	6-9
4502	LARGE	8-10	10-12
4503	X-LARGE	11 & UP	12 & UP



WONDERCUP™

Effective in treating the following foot related pains: bursitis, pressure ulcers, achilles tendonitis, and plantar fasciitis.

PROD #	SIZE	MEN'S SIZES	WOMEN'S SIZES
4504	SMALL	2-4	4-5
4505	MEDIUM	5-7	6-9
4506	LARGE	8-10	10-12
4507	X-LARGE	11 & UP	12 & UP

SILIPOS® BRAND PRODUCTS

WONDERSPORT™



Available in 9 sizes, effective in treating metatarsalgia, diabetic feet, morton's neuroma, fat pad atrophy, elongated metatarsals.

<u>PROD #</u>	<u>MEN'S SIZES</u>	<u>PROD #</u>	<u>WOMEN'S SIZES</u>
4508	5-7	4513	4-5
4509	8-9	4514	6-7
4510	10-11	4515	8-9
4511	12-13	4516	10 & UP
4512	14 & UP		

WONDERSOLE™



Full-length medical-grade silicone provides protection, support, and foot comfort. Softer areas under the heel and metatarsal heads enhance comfort. They can be easily trimmed to fit patient and specific footwear, and are washable and extremely durable.

<u>PRODUCT #</u>	<u>U.S. MEN'S SIZES</u>	<u>U.S. WOMEN'S SIZES</u>
4530	7	5-7
4531	8-9	8-9
4532	10	10-11
4533	11-13	-

HEEL SPUR WITH REMOVABLE PLUG



Provides relief for general or specific heel pain. Features shock-absorbing silicone heel cushion with easily removable plug. May be worn with plug as a comfortable heel pad or without plug for heel relief and protection around the heel spur. It fits all shoe sizes and will not slide in shoe.

<u>PROD#</u>	<u>SIZE</u>	<u>MEN'S SIZES</u>	<u>WOMEN'S SIZES</u>
4001	SMALL	-	5-6.5
4002	MEDIUM	5-7.5	7-9.5
4003	LARGE	8-11.5	10-12

SILIPOS® BRAND PRODUCTS

LOW PROFILE HEEL CUSHION



Features extra soft silicone heel cushion with shock pattern to cushion and diffuse weight. Use with the grid down for maximum shock absorption. Use with the grid up for maximum shear absorption and massaging effect.

PROD#	SIZE	MEN'S SIZES	WOMEN'S SIZES
4004	SMALL	-	5-6.5
4005	MEDIUM	5-7.5	7-9.5
4006	LARGE	8-11.5	10-12

3/4 LENGTH INSOLE



Anatomically designed to support longitudinal and transverse arches. Features a deeper heel cup that provides superior comfort and control. Also has soft, dual-density silicone gel circles providing additional shock absorption and pressure relief to sensitive areas.

PROD#	SIZE	MEN'S SIZES	WOMEN'S SIZES
4010	SMALL	-	5-6.5
4011	MEDIUM	5-7.5	7-9.5
4012	LARGE	8-11.5	10-12

FULL LENGTH INSOLE



Ideal as a replacement insole in shoes with removable insoles. Provides relief for heel, arch, and forefoot problems. Dual formulation absorbs shock under high pressure areas while supporting surrounding tissue. Features multiple vents to enhance air circulation.

PROD#	SIZE	MEN'S SIZES	WOMEN'S SIZES
4013	SMALL	5-6.5	6.5-8.5
4014	MEDIUM	7-9.5	9-11
4015	LARGE	10-12	11.5-12
4016	X-LARGE	12.5	-

SILIPOS® BRAND PRODUCTS

FULL LENGTH INSOLE WITH TOE CREST



Thin design ideal for use in shoes with removable insoles or minimal room. Thin, soft silicone with dual density dots for full-foot comfort and protection. Broad, soft metatarsal raise for relief of forefoot discomfort. Designed to reduce pressure while supporting mid- and forefoot areas. Soft, full toe crest enhances forefoot function.

PROD#	SIZE	MEN'S SIZES	WOMEN'S SIZES
4017	SMALL	-	5-6
4018	MEDIUM	5-6	6.5-7.5
4019	LARGE	6.5-7.5	8-9
4033	X-LARGE	8-9	10-10.5
4034	XX-LARGE	9.5-10	11-11.5



LOW PROFILE FABRIC-COVERED INSOLE

Use for general shock absorption and pressure redistribution. Features soft, durable, and comfortable antimicrobial fabric cover which allows foot to slide easily into shoe. May be used under other insoles for shock absorption and protection.

PROD#	SIZE	MEN'S SIZES	WOMEN'S SIZES
4030	SMALL	-	4.5-7
4031	MEDIUM	5.5-9.5	7.5-11
4032	LARGE	10-12.5	11.5-12



GEL METATARSAL PAD WITH TOE SPREADER™

Anatomically designed to provide comfort and protection of the metatarsal heads. Helps reduce pain and burning sensations under metatarsal heads. Ideal for conditions such as Fat Pad Atrophy, Metatarsalgia, Morton's Neuroma, etc.

PROD#	SIZE	MEN'S SIZES	WOMEN'S SIZES
4020	S/M	5.5-6	4.5-7.5
4021	L/XL	6.5-11.5	8-11

SILIPOS® BRAND PRODUCTS



SILIPOS® DIGITAL™ CARE + KIT

This digital care kit is a convenient way to store, display, and dispense your digital care products. With a special compartment for each product, the care kit allows you to easily select the best solution for the condition you are treating. The kit is discounted and a great value.

PRODUCT # 30365

CONTAINS	SIZE	QTY/ KIT	USES		PROD #
			LOWER EXTREMITY	UPPER EXTREMITY	
DIGITAL PAD	MINI	4			10055
	S/M	4	CORNS	NAIL BED INJURIES	10065
	L/ XL	4			10075
DIGITAL CAP	S/M	6	HAMMERTOES	FINGER TIP INJURIES	10205
MESH TUBING	~3/4" X 6" ~2 cm X 15 cm	2	SORE TOES	SORE FINGERS	10855
GEL TOE SPREADER	MEDIUM	4	STRAIGHTENING/ ALIGNING TOES	STRAIGHTENING/ ALIGNING TOES	11515
	LARGE	4			11525
GEL TOE SEPARATOR	SMALL	6	IRRITATION/ ABRASION	IRRITATION/ ABRASION	11605
	LARGE	6	BETWEEN TOES	BETWEEN TOES	11615

CARPAL GEL PADS™



Same as the pad in the carpal gel sleeve but with self adhesive backing. Can be incorporated into custom made or off-the-shelf hand/ wrist braces.

PRODUCT #	SIZE	SOLD
14205	2" X 6" (5 cm X 15.24 cm)	3 PADS/ PKG

SILIPOS® BRAND PRODUCTS

CARPAL GEL SLEEVE™



A key-hole shaped gel pad molded to the sleeve that fits the contour of the palm. Gel pad designed to cover carpal tunnel scar for placement directly over sutured area. Used to reduce and prevent scar tissue. Insulates against shock and vibration. Helps reduce discomfort caused by hypersensitivity. Can be worn under gloves and/or splints.

PRODUCT #	HAND/ SIZE	SOLD
14115	LEFT/ SMALL	1 PKG
14125	RIGHT/ SMALL	1 PKG
14135	LEFT/ LARGE	1 PKG
14145	RIGHT/ LARGE	1 PKG



HALF COATED BODY SLEEVING™

Nylon/ Spandex Knit Tubing coated on half its circumference with Silipos® Polymer Gel. Can be cut to desired length and will not fray. Gel coating reduces abrasion by moving with the skin. May be used for the following applications: provide compression, protect ulnar plating, reduce scar tissue and increase range of motion as an ulnar nerve sleeve, on donor sites, and for cold therapy.

PRODUCT #	APPROXIMATE SIZE
10895	2 1/4" X 10 1/2' (5.7cm X 3.2m)
10905	2 3/8" X 10 1/2' (6cm X 3.2m)
10915	2 3/4" X 10 1/2' (7cm X 3.2m)
10925	3 1/4" X 10 1/2' (8.3cm X 3.2m)
10935	4" X 10 1/2' (10cm X 3.2m)
10945	5 3/4" X 10 1/2' (14.6cm X 3.2m)

SILIPOS® BRAND PRODUCTS



MESH TUBING™

Elastic tube coated on the inside with gel. Can be worn over ankle and shin or forearm. Also, can be worn under orthosis to reduce friction. Offers gentle compression. Moisturizes dry skin and protects from abrasion.

PRODUCT #	SIZE	SOLD
10845	3" X 10" (7.6 cm X 25.4 cm)	ONE TUBE



GEL-E-ROLL™

An elastic bandage coated on one side with a layer of Silipos® Polymer Gel. Conforms very well to the irregular shapes of the body. Moves with the skin, reducing irritation. May be used on the following applications: reduce scar tissue, for cold therapy, provide compression, residual limb shrinker, neck wrap to treat scar tissue, help aid muscle strains or arthritic pain, moisturize dry skin, in burn care, and to control edema.

PRODUCT #	SIZE	SOLD
10355	4" x 48" (10 cm X 120 cm)	1 ROLL/ PKG



GEL DOTS

1/8" (3mm) x 1" (3cm) diameter dots of gel. A self-adhesive backing allows an individual dot to be adhered to most materials. May be used to reduce abrasion, friction, and pressure points in a splint, cast, brace liner, or footwear. Can be fabricated with insoles under individual metatarsal heads. Use to protect the ulnar styloid or as individual knuckle cushions in a/p splints or circumferential casts.

PRODUCT #	SIZE/ SOLD
10255	FLAT - 15 DOTS/ SHEET

SILIPOS® BRAND PRODUCTS



GEL SQUARES™

3mm thick square of Silipos® Polymer Gel, available with or without self-adhesive backing. A versatile solution to problems caused by pressure and friction. Can be cut with scissors to the required shape and size. Mineral oil in the gel is used to moisturize dry skin. Can be used as an interface between the skin and braces.

PRODUCT #	SIZE	SOLD
15505	SELF-ADHESIVE BACKING 4" X 4" (10 cm X 10 cm)	2 PC./ PKG



PRESSURE RELIEF PADDING™

2mm thick gel padding with a 2-way stretch top cover and adhesive backing. Used to reduce shock and vibration. Can be used as a replacement for moleskin or orthopedic felt in total contact casts.

PRODUCT #	SIZE	SOLD
11205	4" X 36"	1 ROLL/ PKG
11215	8" X 36"	1 ROLL/ PKG



BODY DISCS

Multi purpose discs of Silipos® Polymer Gel. Similar applications as the Gel Squares™. Available in two diameters, 2.5" and 4" (6.3cm and 10cm).

PRODUCT #	SIZE	SOLD
10305	2.5" (6.3 cm DIAMETER)	2/ PKG
10315	4" (10 cm DIAMETER)	2/ PKG


SILIPOS® BRAND PRODUCTS

SOFT SHEAR GEL SHEETING™



3mm thick sheet of Silipos® Polymer Gel, available either lined both sides with fabric or just one side. Fabric cover chosen for stretch qualities. Fabric increases durability of gel, and allows sheet to be glued. Comes with 1 sheet per package.

PRODUCT #	SIZE	FABRIC
15605	8" x 8" (20 cm X 20 cm)	ONE SIDE
15615	8" X 8" (20 cm X 20 cm)	TWO SIDES
15705	12" X 16" (30.5 cm X 40.5 cm)	ONE SIDE
15715	12" X 16" (30.5 cm X 40.5 cm)	TWO SIDES
15805	16" X 20" (40.5 cm X 51 cm)	ONE SIDE
15815	16" X 20" (40.5 cm X 51 cm)	TWO SIDES



WONDERFLEX™ SILICONE SHEETING

Unlined: Recommended to line orthotic appliances.

Lined: Soft durometer cushioning for top cover of orthotics. Ideal for diabetic, arthritic and the insensitive foot. Resists compression and absorbs shock and shear force. Medical grade silicone gel sheeting. May be cut to desired size.

PROD #	SIZE	DESCRIPTION	THICKNESS
4722	12" X 16" (30.5 cm X 40.5 cm)	UNLINED	2MM
4723	12" X 16" (30.5 cm X 40.5 cm)	UNLINED	3MM
4733	12" X 16" (30.5 cm X 40.5 cm)	LINED	3MM
4735	12" X 16" (30.5 cm X 40.5 cm)	LINED	5MM

SILIPOS® BRAND PRODUCTS



SILOSHEATH™ NYLON SHEATHS

Sizing Chart Below For Original, Double Cushion, And Extra Life Silosheath Nylon Sheaths.

SIZE	DESCRIPTION	APPROX. GEL LENGTH	APPROX. LENGTH	APPROX. DISTAL CIRC.*	APPROX. PROXIMA
1	CHILD	8" 20 cm	10" - 13" 25 cm - 33 cm	6.5" - 9.5" 17 cm - 24 cm	6.5" - 10.5" 17 cm - 25 cm
2	SHORT/ NARROW	10" 25 cm	15" - 17" 38 cm - 43 cm	7" - 11" 18 cm - 28 cm	7" - 12.5" 18 cm - 32 cm
3	SHORT/ REGULAR	10" 25 cm	15" - 17" 38 cm - 43 cm	8" - 13" 20 cm - 34 cm	8" - 15" 20 cm - 38 cm
4	MEDIUM/ NARROW	13" 33 cm	17" - 19" 43 cm - 48 cm	8" - 13" 20 cm - 34 cm	8" - 15" 20 cm - 38 cm
5	MEDIUM/ REGULAR	10" 25 cm	17" - 19" 43 cm - 48 cm	9" - 15" 23 cm - 37 cm	9" - 17" 23 cm - 43 cm
6	LARGE/ NARROW	13" 33 cm	18" - 21" 46 cm - 53 cm	9" - 15" 23 cm - 37 cm	9" - 17" 23 cm - 43 cm
7	LARGE/ REGULAR	13" 33 cm	18" - 21" 46 cm - 53 cm	11" - 17" 28 cm - 42 cm	11" - 19" 28 cm - 49 cm
8	LARGE/ WIDE	13" 33 cm	18" - 21" 46 cm - 53 cm	12" - 20" 31 cm - 50 cm	12" - 21" 31 cm - 56 cm

*Please note that when measuring for extra life products, distal and proximal circumferences are approximately 10% less than the above Silosheath™ measurements.



SILOSHEATH™ ORIGINAL

2 ply nylon sheath with 2mm gel.

PRODUCT #	SIZE
13315	1
13325	2
13335	3
13345	4
13355	5
13365	6
13375	7
13385	8

SILIPOS® BRAND PRODUCTS



SILOSHEATH™ DOUBLE CUSHION

2 ply nylon sheath with 2mm gel and extra gel at distal end.

PRODUCT #	SIZE
13415	1
13425	2
13435	3
13445	4
13455	5
13465	6
13475	7
13485	8



SILOSHEATH™ EXTRA LIFE

2 ply nylon sheath with 2mm gel and increased durability.

PRODUCT #	SIZE
13815	1
13825	2
13835	3
13845	4
13855	5
13865	6
13875	7
13885	8

SILIPOS® BRAND PRODUCTS



SOCKET GEL LINERS

SINGLE SOCKET GEL LINER™ - 3 ply gel liner with 3mm gel against the skin, a tapered shape, and Coolmax® exterior material.

PROD #	APPROX. SIZE	APPROX. GEL LENGTH	LENGTH	TOE/ TOP CIRC.
12725	SHORT/ NARROW	10" 25 cm	12" 3.5 cm	9" - 12"/ 15" - 18" 23 cm - 30.5 cm/ 38 cm - 45 cm
12755	SHORT/ REGULAR	10" 25 cm	12" 30.5 cm	12" - 14"/ 15" - 18" 30.5 cm - 35.5 cm/ 38 cm - 45 cm
12765	SHORT/ REGULAR	12" 30.5 cm	16" 40.5 cm	12" - 14"/ 15" - 18" 30.5 cm - 35.5 cm/ 38 cm - 45 cm



DOUBLE SOCKET GEL LINER™ - 5 ply gel liner with 3mm gel and two layers of Coolmax®. Features a tapered shape and has maximum durability.

PROD #	APPROX. SIZE	APPROX. GEL LENGTH	LENGTH	TOE/ TOP CIRC.
12625	SHORT/ NARROW	10" 25 cm	12" 3.5 cm	9" - 12"/ 15" - 18" 23 cm - 30.5 cm/ 38 cm - 45 cm
12655	SHORT/ REGULAR	10" 25 cm	12" 30.5 cm	12" - 14"/ 15" - 18" 30.5 cm - 35.5 cm/ 38 cm - 45 cm
12665	SHORT/ REGULAR	12" 30.5 cm	16" 40.5 cm	12" - 14"/ 15" - 18" 30.5 cm - 35.5 cm/ 38 cm - 45 cm



ULTRA CUSHION GEL LINER™

3 or 5 ply gel liner with 3mm gel and two layers of cotton material. Features a nylon reinforced toe for added durability. This gel liner has the best combination of protection, quality, and value.

PROD #	PROD#	SIZE	APPROX. GEL LENGTH	APPROX. LENGTH OF LINER	APPROX. DISTAL CIRC.	APPROX. CIRC.
12555	12515	I SHORT/ REG.	10" 25 cm	12.5" 32 cm	7" - 11" 18 cm - 28 cm	8" - 13" 20 cm - 33 cm
12565	12525	II MEDIUM / REG.	12" 30.5 cm	13.5" 34 cm	18" - 14" 20 cm - 35.5 cm	10" - 15" 25 cm - 38 cm
12575	12535	III MEDIUM / WIDE	12" 30.5 cm	14.5" 37 cm	10" - 16" 25 cm - 40.5 cm	12" - 17" 30.5 cm - 43 cm

SILIPOS® BRAND PRODUCTS



SYME'S GEL LINER

5 ply gel liner with 3mm vitamin enriched gel, featuring Coolmax® material. Provides optimum conformability and prevents abrasion and skin breakdown. One size fits most.

PRODUCT #	APPROX. LENGTH	APPROX. GEL LENGTH	APPROX. DISTAL CIRC.	APPROX. PROXIMAL CIRC.
1302	16.5" 42 cm	15" 38 cm	10" - 14.5" 25 cm - 37 cm	15" - 20" 38 cm - 51 cm



PARTIAL FOOT SOCK

For use with transmetatarsal amputees. Features elastic conforming Coolmax® material and may be used with foot orthosis or shoes with fillers. Gel covers amputation site.

PRODUCT #	SIZE
10805	REGULAR
10815	X-LONG



MOLDING DISTAL ROLL-ON™

Roll on to distal limb to enhance socket liner fit and combat skin problems. Cut to shape for extra protection anywhere under an Explorer Liner™. Self-adhesive to other Silipos® gel products. Cut into a band of any desired width for extra comfort and protection under the proximal end of a suspension sleeve. Customize the level of protection by using in conjunction with Silosheaths™. Wash with mild soap and water before use.

PRODUCT #	SIZE	APPROX. GEL LENGTH	APPROX. OUTSIDE CIRC.	APPROX. GEL THICKNESS
1495	NARROW	6" / 15 cm	7 1/4" / 15.5 cm	2MM
1496	MEDIUM	6" / 15 cm	8 1/2" / 21.5 cm	2MM
1497	WIDE	6" / 15 cm	9 1/4" / 23.5 cm	2MM