

UNIVERSAL GEL STRAP™

- Surgical-grade elastic sleeve with a gel pad that covers all five metatarsal heads
- Helps reduce pain and discomfort
- Serves as a soft tissue supplement to fat pad atrophy
- Pad is available uncovered for maximum benefit of the mineral oil based gel
- Pad is also available covered to provide an extended life when callusing is not a problem
- By cutting the stitch, can also be used on the instep, wrist or hand

PRODUCT #	SIZE	SOLD
10155	S/M (GEL UNCOVERED)	1/PKG
10165	L/XL (GEL UNCOVERED)	1/PKG
10175	S/M (GEL COVERED)	1/PKG
10185	L/XL (GEL COVERED)	1/PKG

