

## AFO PADS

**ORTHOPEDIC FELT STITCHED TO ORTHOPEDIC VINYL.**



PRODUCT #	FITS STRAP WIDTH	COLOR	DESCRIPTION	PACKED
WBC6124S	5/8"	WHITE	AFO PAD	EACH
WBC6124M	1"	WHITE	AFO PAD	EACH
WBC6124L	1 1/2"	WHITE	AFO PAD	EACH
WBC6124XL	2"	WHITE	AFO PAD	EACH
WBC6124SB	5/8"	BEIGE	AFO PAD	EACH
WBC6124MB	1"	BEIGE	AFO PAD	EACH
WBC6124LB	1 1/2"	BEIGE	AFO PAD	EACH
WBC6124XLB	2"	BEIGE	AFO PAD	EACH
WBC6124SL	5/8"	BLACK	AFO PAD	EACH
WBC6124ML	1"	BLACK	AFO PAD	EACH
WBC6124LL	1 1/2"	BLACK	AFO PAD	EACH
WBC6124XLL	2"	BLACK	AFO PAD	EACH

## SLENDER AFO PADS



PRODUCT #	WIDTH x HEIGHT	COLOR	DESCRIPTION	PACKED
AFO1W	1 3/4" x 2" FELT	WHITE	1" SLENDER AFO PAD	EACH
AFO15W	1 7/8" x 2 5/8" FELT	WHITE	1 1/2" SLENDER AFO PAD	EACH
AFO2W	2 1/2" x 3 3/8" FELT	WHITE	2" SLENDER AFO PAD	EACH
AFO1B	1 3/4" x 2" FELT	BEIGE	1" SLENDER AFO PAD	EACH
AFO15B	1 7/8" x 2 5/8" FELT	BEIGE	1 1/2" SLENDER AFO PAD	EACH
AFO2B	2 1/2" x 3 3/8" FELT	BEIGE	2" SLENDER AFO PAD	EACH
AFO1L	1 3/4" x 2" FELT	BLACK	1" SLENDER AFO PAD	EACH
AFO15L	1 7/8" x 2 5/8" FELT	BLACK	1 1/2" SLENDER AFO PAD	EACH
AFO2L	2 1/2" x 3 3/8" FELT	BLACK	2" SLENDER AFO PAD	EACH

**FOR MATCHING STRAPS,  
SEE PAGES  
64 AND 67**

# AFOs



## AFO PAD AND STRAP KIT

OUR AFO PAD COMPLETE WITH AN ULTRASONICALLY WELDED HOOK AND LOOP STRAP, METAL LINK AND CHAFE FOR EASY ASSEMBLY.

PRODUCT #	WIDTH	COLOR	DESCRIPTION	PACKED
WBC10WD	1"	WHITE	PAD & STRAP KIT	5/PKG
WBC15WD	1½"	WHITE	PAD & STRAP KIT	5/PKG
WBC20WD	2"	WHITE	PAD & STRAP KIT	5/PKG
(AFO pad is white)				
WBC10WDB	1"	BEIGE	PAD & STRAP KIT	5/PKG
WBC15WDB	1½"	BEIGE	PAD & STRAP KIT	5/PKG
WBC20WDB	2"	BEIGE	PAD & STRAP KIT	5/PKG
(AFO pad is white)				
WBC10WDL	1"	BLACK	PAD & STRAP KIT	5/PKG
WBC15WDL	1½"	BLACK	PAD & STRAP KIT	5/PKG
WBC20WDL	2"	BLACK	PAD & STRAP KIT	5/PKG
(AFO pad is white)				



## AFO REINFORCED STRAP KIT

OUR AFO PAD COMPLETE WITH AN ULTRASONICALLY WELDED STRAP WITH HOOK END REINFORCED WITH A PLASTIC WASHER. EACH STRAP ALSO COMES WITH A METAL LINK AND PLASTIC CHAFE FOR EASY ASSEMBLY.

PRODUCT #	WIDTH	COLOR	DESCRIPTION	PACKED
6101	1"	WHITE	15" STRAP KIT	5/PKG
6115	1½"	WHITE	15" STRAP KIT	5/PKG
6120	2"	WHITE	15" STRAP KIT	5/PKG
(AFO pad is white)				
6201	1"	BEIGE	15" STRAP KIT	5/PKG
6215	1½"	BEIGE	15" STRAP KIT	5/PKG
6220	2"	BEIGE	15" STRAP KIT	5/PKG
(AFO pad is beige*)				
6301	1"	BLACK	15" STRAP KIT	5/PKG
6315	1½"	BLACK	15" STRAP KIT	5/PKG
6320	2"	BLACK	15" STRAP KIT	5/PKG
(AFO pad is black*)				

\* white pads are available if specified.