

# STRAPS



## WASHER REINFORCED HOOK & LOOP STRAPS

With a plastic washer ultrasonically welded into hook end and completely encased in the material, these straps are great for AFOs and KAFO's.

PRODUCT #	WIDTH/LENGTH	COLOR	PACKAGE
112WWASH	1" x 12"	WHITE STRAPS	25/PKG
116WWASH	1" x 16"	WHITE STRAPS	25/PKG
120WWASH	1" x 20"	WHITE STRAPS	25/PKG
1512WWASH	1 1/2" x 12"	WHITE STRAPS	25/PKG
1516WWASH	1 1/2" x 16"	WHITE STRAPS	25/PKG
1520WWASH	1 1/2" x 20"	WHITE STRAPS	25/PKG
212WWASH	2" x 12"	WHITE STRAPS	25/PKG
216WWASH	2" x 16"	WHITE STRAPS	25/PKG
220WWASH	2" x 20"	WHITE STRAPS	25/PKG
112BWASH	1" x 12"	BEIGE STRAPS	25/PKG
116BWASH	1" x 16"	BEIGE STRAPS	25/PKG
120BWASH	1" x 20"	BEIGE STRAPS	25/PKG
1512BWASH	1 1/2" x 12"	BEIGE STRAPS	25/PKG
1516BWASH	1 1/2" x 16"	BEIGE STRAPS	25/PKG
1520BWASH	1 1/2" x 20"	BEIGE STRAPS	25/PKG
212BWASH	2" x 12"	BEIGE STRAPS	25/PKG
216BWASH	2" x 16"	BEIGE STRAPS	25/PKG
220BWASH	2" x 20"	BEIGE STRAPS	25/PKG
112LWASH	1" x 12"	BLACK STRAPS	25/PKG
116LWASH	1" x 16"	BLACK STRAPS	25/PKG
120LWASH	1" x 20"	BLACK STRAPS	25/PKG
1512LWASH	1 1/2" x 12"	BLACK STRAPS	25/PKG
1516LWASH	1 1/2" x 16"	BLACK STRAPS	25/PKG
1520LWASH	1 1/2" x 20"	BLACK STRAPS	25/PKG
212LWASH	2" x 12"	BLACK STRAPS	25/PKG
216LWASH	2" x 16"	BLACK STRAPS	25/PKG
220LWASH	2" x 20"	BLACK STRAPS	25/PKG

# STRAPS



## SELF ADHESIVE HOOK & LOOP STRAPS

1" X 4" PRESSURE SENSITIVE HOOK ULTRASONICALLY WELDED TO SEW QUALITY LOOP IN CONVENIENT READY TO USE PACKAGES. JUST PEEL OFF LINER ON THE HOOK AND PLACE WHERE NEEDED, THEN TRIM THE LOOP TO DESIRED LENGTH.

PRODUCT #	WIDTH & LENGTH	COLOR	PACKED
114WSPS	1" x 14"	WHITE	25/PKG
214WSPS	2" x 14"	WHITE	25/PKG
114BSPS	1" x 14"	BEIGE	25/PKG
214BSPS	2" x 14"	BEIGE	25/PKG
114LSPS	1" x 14"	BLACK	25/PKG
214LSPS	2" x 14"	BLACK	25/PKG
114RBSPS	1" x 14"	ROYAL BLUE	25/PKG
214RBSPS	2" x 14"	ROYAL BLUE	25/PKG



## SELF ADHESIVE HOOK & PLUSH PILE STRAPS


1" X 3" PRESSURE SENSITIVE HOOK ULTRASONICALLY WELDED TO FOAM PADDED LOOP IN CONVENIENT READY TO USE PACKAGES. JUST PEEL OFF LINER ON THE HOOK AND PLACE WHERE NEEDED, THEN TRIM THE LOOP TO DESIRED LENGTH. EXTREMELY DURABLE, SOFT WITH NO STRETCH.

PRODUCT #	WIDTH & LENGTH	COLOR	PACKED
PP116W	1" x 16"	WHITE	25/PKG
PP216W	2" x 16"	WHITE	25/PKG
PP116B	1" x 16"	BEIGE	25/PKG
PP216B	2" x 16"	BEIGE	25/PKG


**LOOKING FOR AFO PADS?**

**CHECK OUT WBC'S LARGE SELECTION OF AFO PADS**

**GO TO PAGE 73**



# STRAPS



## PLUSH PILE D-RING STRAPS

PADDED LOOP MATERIAL WITH A STRONG 2" D-RING.  
SOFT FOAM REDUCES SKIN IRRITATION.

PRODUCT#	WIDTH & LENGTH	COLOR	PACKED
DPP214W	2" x 14"	WHITE	25/PKG
DPP214B	2" x 14"	BEIGE	25/PKG




## HOOK & LOOP D-RING STRAPS

SIMPLY THREAD THE HOOK END THROUGH THE D-RING AND  
ATTACH TO LOOP MATERIAL.

PRODUCT #	WIDTH & LENGTH	COLOR	PACKED
D112W	1" x 12"	WHITE	25/PKG
D116W	1" x 16"	WHITE	25/PKG
D218W	2" x 18"	WHITE	25/PKG

***DID YOU KNOW?***  
***FULL CASES***  
***OF PLASTAZOTE SHIP***  
***FREE OF DELIVERY CHARGE!***



# STRAPS

## APLIX® SELF GRIPPING FASTENER STRAPS

### WASHER REINFORCED STRAPS


MADE WITH APLIX® BRAND UNNAPPED TEXTURIZED LOOP.



PRODUCT #	WIDTH & LENGTH	COLOR	PACKED
120WWASHAP	1" x 20"	WHITE	25/PKG
1520WWASHAP	1½" x 20"	WHITE	25/PKG
220WWASHAP	2" x 20"	WHITE	25/PKG
120BWASHAP	1" x 20"	BEIGE	25/PKG
1520BWASHAP	1½" x 20"	BEIGE	25/PKG
220BWASHAP	2" x 20"	BEIGE	25/PKG
120LWASHAP	1" x 20"	BLACK	25/PKG
1520LWASHAP	1½" x 20"	BLACK	25/PKG
220LWASHAP	2" x 20"	BLACK	25/PKG

### 50/50 APLIX® BRAND STRAPS

MADE WITH HALF HOOK AND HALF APLIX® BRAND UNNAPPED TEXTURIZED LOOP.



PRODUCT #	WIDTH & LENGTH	COLOR	PACKED
H112WSAP	1" x 12"	WHITE	25/PKG
H116WSAP	1" x 16"	WHITE	25/PKG
H1512WSAP	1½" x 12"	WHITE	25/PKG
H1516WSAP	1½" x 16"	WHITE	25/PKG
H214WSAP	2" x 14"	WHITE	25/PKG
H218WSAP	2" x 18"	WHITE	25/PKG
H112BSAP	1" x 12"	BEIGE	25/PKG
H116BSAP	1" x 16"	BEIGE	25/PKG
H1512BSAP	1½" x 12"	BEIGE	25/PKG
H1516BSAP	1½" x 16"	BEIGE	25/PKG
H214BSAP	2" x 14"	BEIGE	25/PKG
H218BSAP	2" x 18"	BEIGE	25/PKG
H112LSAP	1" x 12"	BLACK	25/PKG
H116LSAP	1" x 16"	BLACK	25/PKG
H1512LSAP	1½" x 12"	BLACK	25/PKG
H1516LSAP	1½" x 16"	BLACK	25/PKG
H214LSAP	2" x 14"	BLACK	25/PKG
H218LSAP	2" x 18"	BLACK	25/PKG

# STRAPS




## DELUXE "C" FOLD ADJUSTABLE STRAP KITS

OUR "PATIENT-FRIENDLY" C-FOLD WITH ROUND HOOK SEWN ON ONE END. EACH STRAP COMES WITH A METAL LINK AND CHAFE. WE THEN WELD THE OTHER END WITH 3 CUTTING LINES FOR EASY SIZING ADJUSTMENTS WITHOUT FRAYING.

PRODUCT #	WIDTH & LENGTH	COLOR/DESCRIPTION	PACKED
DCFS1010W	1" x 10"	WHITE STRAP KIT	5/PKG
DCFS1514W	1½" x 14"	WHITE STRAP KIT	5/PKG
DCFS1518W	1½" x 18"	WHITE STRAP KIT	5/PKG
DCFS2018W	2" x 18"	WHITE STRAP KIT	5/PKG
DCFS1010B	1" x 10"	BEIGE STRAP KIT	5/PKG
DCFS1514B	1½" x 14"	BEIGE STRAP KIT	5/PKG
DCFS1518B	1½" x 18"	BEIGE STRAP KIT	5/PKG
DCFS2018B	2" x 18"	BEIGE STRAP KIT	5/PKG

**NOTE:** EACH PACKAGE CONTAINS FIVE (5) ADJUSTABLE STRAPS PER BAG. EACH STRAP HAS A METAL LINK WITH PLASTIC CHAFE.



## DOUBLE FACED KNIT LOOP STRAPS

SOFT KNIT LOOP, HOOK COMPATIBLE ON BOTH SIDES. MAKES A GREAT STRAP, PROVIDING STRENGTH AND COMFORT. WE SEW ON INDUSTRIAL HOOK FOR A GREAT HOLD.

PRODUCT #	WIDTH & LENGTH	COLOR	PACKED
DF1518SW	1½" x 18"	WHITE	25/PKG
DF218SW	2" x 18"	WHITE	25/PKG
DF1518SB	1½" x 18"	BEIGE	25/PKG
DF218SB	2" x 18"	BEIGE	25/PKG
DF1518SL	1½" x 18"	BLACK	25/PKG
DF218SL	2" x 18"	BLACK	25/PKG


# STRAPS



## ZIG-ZAG ANGLED STRAPS

LOOP ON BOTH SIDES WITH HOOK ON HALF OF ONE SIDE. STITCHED CLOSED WITH ZIG-ZAG PATTERN FOR CUT-OFF STRIPS AT ANY POINT.

PRODUCT #	WIDTH/LENGTH	COLOR	DESCRIPTION	QUANTITY
ZZ118W	1" x 18"	WHITE	ANGLED STRAP	10/PKG
ZZ1518W	1½" x 18"	WHITE	ANGLED STRAP	10/PKG
ZZ118B	1" x 18"	BEIGE	ANGLED STRAP	10/PKG
ZZ1518B	1½" x 18"	BEIGE	ANGLED STRAP	10/PKG
ZZ118L	1" x 18"	BLACK	ANGLED STRAP	10/PKG
ZZ1518L	1½" x 18"	BLACK	ANGLED STRAP	10/PKG




## DACRON® BACKED HOOK & LOOP STRAPS

HIGH QUALITY WHITE DACRON STITCHED WITH HOOK & LOOP.  
***GREAT FOR BODY JACKETS.***

PRODUCT #	WIDTH/LENGTH	COLOR	PACKAGED
1510WDS	1½" x 10"	WHITE	10/PKG
1512WDS	1½" x 12"	WHITE	10/PKG
1516WDS	1½" x 16"	WHITE	10/PKG
1518WDS	1½" x 18"	WHITE	10/PKG
1520WDS	1½" x 20"	WHITE	10/PKG
1522WDS	1½" x 22"	WHITE	10/PKG
210WDS	2" x 10"	WHITE	10/PKG
212WDS	2" x 12"	WHITE	10/PKG
216WDS	2" x 16"	WHITE	10/PKG
218WDS	2" x 18"	WHITE	10/PKG
220WDS	2" x 20"	WHITE	10/PKG
222WDS	2" x 22"	WHITE	10/PKG



# STRAPS



## WBC BRAND 25/75 STRAPS

25% HOOK AND BALANCE LOOP WELDED TOGETHER.

PRODUCT #	WIDTH & LENGTH	COLOR	PACKED
120WFS	1" x 20"	WHITE	25/PKG
124WFS	1" x 24"	WHITE	25/PKG
1520WFS	1½" x 20"	WHITE	25/PKG
1524WFS	1½" x 24"	WHITE	25/PKG
220WFS	2" x 20"	WHITE	25/PKG
224WFS	2" x 24"	WHITE	25/PKG
120BFS	1" x 20"	BEIGE	25/PKG
124BFS	1" x 24"	BEIGE	25/PKG
1520BFS	1½" x 20"	BEIGE	25/PKG
1524BFS	1½" x 24"	BEIGE	25/PKG
220BFS	2" x 20"	BEIGE	25/PKG
224BFS	2" x 24"	BEIGE	25/PKG
120LFS	1" x 20"	BLACK	25/PKG
124LFS	1" x 24"	BLACK	25/PKG
1520LFS	1½" x 20"	BLACK	25/PKG
1524LFS	1½" x 24"	BLACK	25/PKG
220LFS	2" x 20"	BLACK	25/PKG
224LFS	2" x 24"	BLACK	25/PKG

## WBC BRAND 50/50 STRAPS

HALF HOOK AND HALF LOOP WELDED TOGETHER WITH ROUND END ON HOOK SIDE.




PRODUCT #	WIDTH & LENGTH	COLOR	PACKED
H112WS	1" x 12"	WHITE	25/PKG
H116WS	1" x 16"	WHITE	25/PKG
H120WS	1" x 20"	WHITE	25/PKG
H1512WS	1½" x 12"	WHITE	25/PKG
H1516WS	1½" x 16"	WHITE	25/PKG
H1520WS	1½" x 20"	WHITE	25/PKG
H214WS	2" x 14"	WHITE	25/PKG
H218WS	2" x 18"	WHITE	25/PKG
H222WS	2" x 22"	WHITE	25/PKG
H224WS	2" x 24"	WHITE	25/PKG

# STRAPS



## WBC BRAND 50/50 STRAPS (continued)

PRODUCT #	WIDTH & LENGTH	COLOR	PACKED
H112BS	1" x 12"	BEIGE	25/PKG
H116BS	1" x 16"	BEIGE	25/PKG
H120BS	1" x 20"	BEIGE	25/PKG
H1512BS	1½" x 12"	BEIGE	25/PKG
H1516BS	1½" x 16"	BEIGE	25/PKG
H1520BS	1½" x 20"	BEIGE	25/PKG
H214BS	2" x 14"	BEIGE	25/PKG
H218BS	2" x 18"	BEIGE	25/PKG
H222BS	2" x 22"	BEIGE	25/PKG
H224BS	2" x 24"	BEIGE	25/PKG
H112LS	1" x 12"	BLACK	25/PKG
H116LS	1" x 16"	BLACK	25/PKG
H120LS	1" x 20"	BLACK	25/PKG
H1512LS	1½" x 12"	BLACK	25/PKG
H1516LS	1½" x 16"	BLACK	25/PKG
H1520LS	1½" x 20"	BLACK	25/PKG
H214LS	2" x 14"	BLACK	25/PKG
H218LS	2" x 18"	BLACK	25/PKG
H222LS	2" x 22"	BLACK	25/PKG
H224LS	2" x 24"	BLACK	25/PKG



## VINYL BACKED HOOK & LOOP STRAPS

SOFT ORTHOPEDIC VINYL-BACKED HOOK & LOOP STRAPS, ELECTRONICALLY SEALED AND BACKED WITH A WATERPROOF, MARINE GRADE VINYL. SOFT, FLEXIBLE AND VERY PATIENT FRIENDLY.

PRODUCT #	SIZE	COLOR	PACKED
LLS118L	1" x 18"	BLACK	10/PACK
LLS118W	1" x 18"	WHITE	10/PACK
LLS118SS	1" x 18"	SANDSTONE	10/PACK
LLS1518L	1½" x 18"	BLACK	10/PACK
LLS1518W	1½" x 18"	WHITE	10/PACK
LLS1518SS	1½" X18"	SANDSTONE	10/PACK
LLS218L	2" x 18"	BLACK	10/PACK
LLS218W	2" x 18"	WHITE	10/PACK
LLS218SS	2" x 18"	SANDSTONE	10/PACK

# STRAPS

## CUSTOM HOOK & LOOP STRAPS

### QUESTIONS & ANSWERS




- Q.** Can I get custom straps to my exact specifications?
- A.** YES, YOU CAN HAVE THE LENGTH AND WIDTH YOU WANT, WELDED OR SEWN, WITH A VARIETY OF HARDWARE, FASTENERS, AND WEBBINGS AVAILABLE.

- Q.** What about colors? Can I get purple, red, blue, or even neon colors?
- A.** YES, WITH OVER 25 DIFFERENT COLORS IN STOCK, YOU'RE SURE TO FIND THE COLOR YOU NEED!



- Q.** I really don't need thousands of straps. What is the minimum order?
- A.** AT WBC INDUSTRIES, THERE ARE NO MINIMUMS ON ORDERS! GIVE US A CALL AND LET US SHOW YOU THE ADVANTAGES OF DEALING WITH A SMALL, FAMILY-OWNED BUSINESS.



- Q.** I'm tired of waiting weeks for straps I need. When can you deliver?
- A.** DELIVERY AT WBC INDUSTRIES IS SAME-DAY SHIPPING ON STOCK STRAPS. CUSTOM ORDERS GENERALLY TAKE THREE TO FIVE WORKING DAYS.
- Q.** I want straps made with top quality materials and rigid quality control standards. Am I going to get all that?
- A.** YOU BET! WBC INDUSTRIES HOOK AND LOOP STRAPS ARE TOP QUALITY. SINCE ALL FABRICATING IS DONE IN-HOUSE, WE HAVE TOTAL CONTROL OVER QUALITY AND COSTS.