

AFO PADS

AFO PADS

Orthopedic felt stitched to orthopedic vinyl.



PRODUCT #	FITS STRAP WIDTH	COLOR	DESCRIPTION	PACKED
WBC6124S	5/8"	WHITE	AFO PAD	EACH
WBC6124M	1"	WHITE	AFO PAD	EACH
WBC6124L	1 1/2"	WHITE	AFO PAD	EACH
WBC6124XL	2"	WHITE	AFO PAD	EACH
WBC6124SB	5/8"	BEIGE	AFO PAD	EACH
WBC6124MB	1"	BEIGE	AFO PAD	EACH
WBC6124LB	1 1/2"	BEIGE	AFO PAD	EACH
WBC6124XLB	2"	BEIGE	AFO PAD	EACH
WBC6124SL	5/8"	BLACK	AFO PAD	EACH
WBC6124ML	1"	BLACK	AFO PAD	EACH
WBC6124LL	1 1/2"	BLACK	AFO PAD	EACH
WBC6124XLL	2"	BLACK	AFO PAD	EACH

SLENDER AFO PADS

Orthopedic felt stitched to orthopedic vinyl. Slit openings are vertical for a slender fit.



PROD #	WIDTH X HEIGHT	COLOR	DESCRIPTION	PACKED
AFO1W	1 3/4" X 2" FELT	WHITE	1" SLENDER AFO PAD	EACH
AFO15W	1 7/8" X 2 5/8" FELT	WHITE	1 1/2" SLENDER AFO PAD	EACH
AFO2W	2 1/2" X 3 3/8" FELT	WHITE	2" SLENDER AFO PAD	EACH
AFO1B	1 3/4" X 2" FELT	BEIGE	1" SLENDER AFO PAD	EACH
AFO15B	1 7/8" X 2 5/8" FELT	BEIGE	1 1/2" SLENDER AFO PAD	EACH
AFO2B	2 1/2" X 3 3/8" FELT	BEIGE	2" SLENDER AFO PAD	EACH
AFO1L	1 3/4" X 2" FELT	BLACK	1" SLENDER AFO PAD	EACH
AFO15L	1 7/8" X 2 5/8" FELT	BLACK	1 1/2" SLENDER AFO PAD	EACH
AFO2L	2 1/2" X 3 3/8" FELT	BLACK	2" SLENDER AFO PAD	EACH

AFO PADS

DOUBLE THICK AFO PADS

Orthopedic felt 1/4" thick stitched to orthopedic vinyl.



PRODUCT #	FITS STRAP		DESCRIPTION	PACKED
	WIDTH	COLOR		
WBC6124S2	5/8"	WHITE	DOUBLE THICK AFO PAD	EACH
WBC6124M2	1"	WHITE	DOUBLE THICK AFO PAD	EACH
WBC6124L2	1 1/2"	WHITE	DOUBLE THICK AFO PAD	EACH
WBC6124XL2	2"	WHITE	DOUBLE THICK AFO PAD	EACH

AFO PAD AND STRAP KIT

Our AFO pad complete with an ultrasonically welded hook & loop strap, metal link, and chafe for easy assembly.



PRODUCT #	WIDTH	COLOR	DESCRIPTION	PACKED
WBC10WD	1"	WHITE	PAD & STRAP KIT	5/ PKG
WBC15WD	1 1/2"	WHITE	PAD & STRAP KIT	5/ PKG
WBC20WD	2"	WHITE	PAD & STRAP KIT	5/ PKG

AFO Pad is White

WBC10WDB	1"	BEIGE	PAD & STRAP KIT	5/ PKG
WBC15WDB	1 1/2"	BEIGE	PAD & STRAP KIT	5/ PKG
WBC20WDB	2"	BEIGE	PAD & STRAP KIT	5/ PKG

AFO Pad is Beige (White Also Available Upon Request)

WBC10WDL	1"	BLACK	PAD & STRAP KIT	5/ PKG
WBC15WDL	1 1/2"	BLACK	PAD & STRAP KIT	5/ PKG
WBC20WDL	2"	BLACK	PAD & STRAP KIT	5/ PKG

AFO Pad is Black (White Also Available Upon Request)

AFO PADS

AFO REINFORCED STRAP KIT



Our AFO pad complete with an ultrasonically welded strap with the hook end reinforced with a strong plastic washer. Each strap also comes with a metal link and plastic chafe for easy assembly.

PRODUCT #	WIDTH	COLOR	DESCRIPTION	PACKED
6101	1"	WHITE	15" STRAP KIT	5/ PKG
6115	1 1/2"	WHITE	15" STRAP KIT	5/ PKG
6120	2"	WHITE	15" STRAP KIT	5/ PKG

AFO Pad is White

6201	1"	BEIGE	15" STRAP KIT	5/ PKG
6215	1 1/2"	BEIGE	15" STRAP KIT	5/ PKG
6220	2"	BEIGE	15" STRAP KIT	5/ PKG

AFO Pad is Beige (White Also Available Upon Request)

6301	1"	BLACK	15" STRAP KIT	5/ PKG
6315	1 1/2"	BLACK	15" STRAP KIT	5/ PKG
6320	2"	BLACK	15" STRAP KIT	5/ PKG

AFO Pad is Black (White Also Available Upon Request)

ADDITIONAL COLORS NOW AVAILABLE FOR STRAP KITS:

Red, Yellow, Kelly Green, Purple, Orange, Royal Blue, Light Gray, Slate Gray, Light Blue, Pink, Brown, Navy Blue, Olive Green, Neon Yellow, Neon Orange, Neon Green, Neon Pink, Burgundy.



AFO TIBIAL PADS WITH COMPLETE CLOSURE KIT

This kit comes complete with AFO pad, strap, chafe, and rivets. The straps are 18" long (half hook and half loop) with rivet hole in loop end. Chafes are available with plastic loops, rectangular metal loops, or metal loops with roller rings.



PROD #	WIDTH	COLOR	DESCRIPTION	PACKAGE
40009	1"	WHITE	PAD, STRAP, RIVETS, CHAFE WITH PLASTIC LOOP	1 SET
40010	1.5"	WHITE	PAD, STRAP, RIVETS, CHAFE WITH PLASTIC LOOP	1 SET
40011	2"	WHITE	PAD, STRAP, RIVETS, CHAFE WITH PLASTIC LOOP	1 SET

Continued →

AFO PADS

AFO TIBIAL PADS WITH COMPLETE CLOSURE KIT, Cont.



PROD #	WIDTH	COLOR	DESCRIPTION	PACKAGE
11009	1"	WHITE	PAD, STRAP, RIVETS, CHAFE WITH RECTANGULAR LOOP	1 SET
11010	1.5"	WHITE	PAD, STRAP, RIVETS, CHAFE WITH RECTANGULAR LOOP	1 SET
11011	2"	WHITE	PAD, STRAP, RIVETS, CHAFE WITH RECTANGULAR LOOP	1 SET
110090	1"	WHITE	PAD, STRAP, RIVETS, CHAFE WITH ROLLER RINGS	1 SET
110100	1.5"	WHITE	PAD, STRAP, RIVETS, CHAFE WITH ROLLER RINGS	1 SET
110110	2"	WHITE	PAD, STRAP, RIVETS, CHAFE WITH ROLLER RINGS	1 SET