

STOCKINETTE



COTTON STOCKINETTE *By Comfort Products*

High quality stockinette made using only 100% natural cotton. Available in 25 yard rolls in a convenient dispenser box.

PRODUCT #	DESCRIPTION	SIZE
SCC1	COTTON STOCKINETTE	1"
SCC15	COTTON STOCKINETTE	1.5"
SCC2	COTTON STOCKINETTE	2"
SCC25	COTTON STOCKINETTE	2.5"
SCC3	COTTON STOCKINETTE	3"
SCC4	COTTON STOCKINETTE	4"
SCC5	COTTON STOCKINETTE	5"
SCC6	COTTON STOCKINETTE	6"
SCC8	COTTON STOCKINETTE	8"
SCC10	COTTON STOCKINETTE	10"
SCC12	COTTON STOCKINETTE	12"



NYLON STOCKINETTE *By Comfort Products*

High quality stockinette made from 100% virgin nylon, which helps achieve a nice full body with a fine gauge rib knit finish. Unlike other stockinettes, this one features greater elasticity and a smoother finish. Comes in 25 yard rolls in a convenient dispenser box.

PRODUCT #	DESCRIPTION	SIZE
SCN1	NYLON STOCKINETTE	1"
SCN2	NYLON STOCKINETTE	2"
SCN25	NYLON STOCKINETTE	2.5"
SCN3	NYLON STOCKINETTE	3"
SCN4	NYLON STOCKINETTE	4"
SCN5	NYLON STOCKINETTE	5"
SCN6	NYLON STOCKINETTE	6"
SCN8	NYLON STOCKINETTE	8"
SCN10	NYLON STOCKINETTE	10"
SCN12	NYLON STOCKINETTE	12"

STOCKINETTE



MULTI-SIZE® NYLON STOCKINETTE *By Comfort Products*

Unique tubular stockinette that has multiple uses. Exceptionally smooth and does not run from either direction. Features a remarkable stretch capacity that allows 4 sizes to accommodate most orthotic and prosthetic requirements. Ideal for final layer laminations. Great choice for vacuum forming for a smooth surface and can be used effectively for flexible polyurethane and similar type spray coatings. Also, can accommodate rapid size changes without folds or creases. Comes in 25 yard rolls.

PROD #	SIZE	FLAT WIDTH	MIN. CIRC.	MAX CIRC.
MSS	SMALL	3"	4"	21"
MSR	REGULAR	4"	6"	25"
MSL	LARGE	5"	7"	31"
MSB	BODY	10"	20"	60"



SPECTRACARB™ STOCKINETTE *By Comfort Products*

Unique tubular stockinette featuring a combination of Spectra® and Carbon. Spectra® fiber is one of the world's strongest and lightest fibers. Most commonly used in lay-ups that alternate two layers of SpectraCarb™ with two layers of Multi-Size® Nylon. Material does not itch and has heat-and-relieve characteristics that allow changes in the laminated device with use of a heat gun.

PROD #	SIZE	DESCRIPTION	ROLL SIZE
CSC2	2"	SPECTRACARB™	5 YARDS
CSC3	3"	SPECTRACARB™	5 YARDS
CSC4	4"	SPECTRACARB™	5 YARDS
CSC5	5"	SPECTRACARB™	5 YARDS
CSC6	6"	SPECTRACARB™	5 YARDS
CSC8	8"	SPECTRACARB™	5 YARDS